



**TRIO.BAS<sup>TM</sup>**  
BIOLOGICAL AIR SAMPLER



# MICROBIAL ENVIRONMENTAL MONITORING 2025

**orum**  
Orum International

# A WIDE RANGE OF SOLUTIONS TO FIT YOUR BEST PRACTICE



CALIBRATION  
CERTIFICATE  
ISSUED WITH EACH  
INSTRUMENT



## YELLOW FAMILY







**TRIO.BAS<sup>TM</sup>**  
BIOLOGICAL AIR SAMPLER

**45**  
YEARS OF AIR  
MONITORING 1976  
2021

**STAINLESS  
STEEL  
FAMILY**



# QUESTIONS AND ANSWERS IS DEDICATED TO THE INTRODUCTION OF MICROBIAL AIR MONITORING (MAM)

## 1 - WHICH IS THE TERMINOLOGY (GLOSSARY) INVOLVED IN MAM?

Active sampling, Air flow l/m, ALCOA, Bacteria, Colony forming unit (Cfu), Contact plate, Contamination, Culture Media, Culture plate, Data integrity, Disinfection, Fungi, Incubation, PQ, IQ, OQ Document, ISO, Yeast, Micron, Micro-organism, Moulds, Passive sampling, Petri dish, SDA agar media, Settle, Sterilization, TSA agar media, USP, Unviable particles, Viable particles

## 2 - WHY IS THE MICROBIAL CONTENT OF THE AIR TESTED?

The air contains micro-organisms that are "positioned" on the airborne dust. Such micro-organisms can contaminate sterile pharmaceutical products, reduce the shelf-life of food, dairy products, infect patients in hospital, contaminate by moulds HVAC of building, destroy old artistic pictures, contaminate theatres, offices, schools, trains, etc.

## 3 - WHAT'S THE DIFFERENCE BETWEEN A PARTICLE COUNTER AND A MICROBIAL AIR SAMPLER?

The particle counter takes a quickly measuring of unviable particles; the microbial air sampler takes a measuring of viable particles. Particle counter cannot distinguish viable and unviable particles and doesn't provide information on the microbiological content of the environment. The particle counter is the answer for an immediate unviable particle result, the microbial air sampler is the answer for a microbial content after incubation. Both methods are therefore used for E.M.

## 4 - WHAT IS THE DEFINITION OF MAM (MICROBIOLOGICAL AIR MONITORING)?

Microbial Air Monitoring (also called Viable Air Monitoring) has the purpose to capture samples of micro-organisms (bacteria, moulds, yeast) suspended in the air on a nutrient media (inside a Petri dish) with subsequent incubation to highlight their growth after 2-3 days for bacteria and 3-6 days for moulds and yeasts.

## 5 - WHAT TYPE OF NUTRIENT AGAR IS USED FOR MAM?

Tryptic Soy Agar (TSA) is typically used for total bacterial count of bacteria; Sabouraud Dextrose Agar (SDA) is used for mould or yeast growth.

The minimum amount of agar in each Petri dish should be 12-15 ml for Contact plate and 25 -30 – 35 ml for standard Petri dish.

## 6 - WHAT IS THE INCUBATION TIME AND TEMPERATURE?

Bacterial samples are typically incubated for 48-72 hours at 30-37°C. Moulds and yeasts are typically incubated for 4-6 days at 26-30°C.

## 7 - WHAT IS THE DIFFERENCE BETWEEN ACTIVE MICROBIAL AIR SAMPLING AND PASSIVE (SETTLE) METHODS?

Active microbial air sampler aspirates a known volume of air that impacts at constant speed on the agar media surface collecting viable particles of different size in microns.

With Passive method (settle) the viable particles fall naturally onto the surface of an open Petri dish and their collection vary based on air flow, humidity, temperature and size of micro-organisms.

The active method reports quantitative results in the format of cfu / 1000 litres of air (=1 cubic meter); the settle method reports the results that fall on open agar surface in 4 hours.

The number of cfu collected by the active method is much larger than with settle method and the quantitative results can be applied to action and alert levels.

## 8 - WHICH IS THE FREQUENCY OF MICROBIAL AIR SAMPLING?

The frequency varies by industry and application.

The pharma companies that produce sterile products have to follow the indication of the Health Authorities and the CCS (Contamination Control System) specific for each producer.

Indications can be found in the GMP Annex 1.

The frequency is daily for high-risk area and bi-weekly or weekly for other sites.

The food, dairy, agro, poultry, beverage companies should make a preliminary investigation in their premises and production sites to find the best interval (daily, weekly, monthly?). Dedicated Application Notes are available on request.

The frequency in hospital varies from high-risk sites (e.g.: type of operating room), high risk patient wards, etc, some indications are reported in the USP 1116 document Microbiological evaluation of Cleanrooms and other controlled environments.

Compounding pharmacies are requested a calibration every 6 months.

In the case of building, the HVAC should be monitored inside and outside to compare the results as number of cfu. The inside number should never be higher than outside.



## **9 - HOW OFTEN AN ACTIVE MICROBIAL AIR SAMPLER SHOULD BE CALIBRATED?**

The time is bounded to the frequency of use, but generally is every 6-12 months.

## **10 - SHOULD THE AIR FLOW RATE OF AN ACTIVE AIR SAMPLER CONTROL IN THE PERIOD BETWEEN THE OLD AND THE NEW CALIBRATION CHECK?**

Yes. If the calibration is out of limit during the final check, all the results of the previous period are not accepted. A check at regular interval is therefore requested using the specific instrument.

## **11 - WHICH IS THE CORRECT VOLUME OF AIR TO BE ASPIRATED BY AN ACTIVE AIR SAMPLER?**

The volume depends from the expected number of micro-organisms present in the environment under control. The number of cfu /plate should be easy to count (not more than about 150 cfu) and therefore it is necessary to establish the volume to be aspirated (e.g.: 300 – 500 – 1000 litres). The obtained result should be then related to 1000 litres of air (1 cubic metre of air).

## **12 - SHOULD THE CFU OBTAINED FINAL RESULT OF AN ACTIVE AIR SAMPLER MODIFIED BY A STATISTICAL TABLE?**

The instruction manual of each sampler reports a statistical table that should be applied to overcome the problem that more than a micro-organism is captured on the same hole of the agar surface.

## **13 - WHY ARE STERILE ASPIRATING HEAD FOR DAILY USE AVAILABLE?**

The aspirating head of the air sampler is typically stainless steel made and must be autoclaved each day with related documentation. The availability of sterile techno-polymer heads (Daily Shift) that are used for a complete daily working shift avoids all the sterilisation iter and eliminated the possible false results.

## **14 - WHICH ACTIVE AIR SAMPLER SHOULD BE USED IN ISOLATOR AND RABS?**

The command unit should be outside from the isolator to save inside space and reduce the risk of microbial contamination. The only part inside is the aspirating head or the funnel model.

## **15 - HOW ARE COMPRESSED AIR / GAS MONITORED FOR POSSIBLE CONTAMINATION IN GRADE A AND B?**

The active air sampler connected to the TRIO.GAS instrument is used according to the guidelines outlined on the ISO 8573-7 document.

## **16 - IS A RE-VALIDATION NECESSARY FOR A NEW ACTIVE AIR SAMPLER-INSTRUMENT INTRODUCTION?**

Yes. A New validation is necessary. A specific protocol is available.

## **17 - WHY TO USE AN ACTIVE AIR SAMPLER WITH AN ASPIRATION OF 200 LITRES PER MINUTE INSTEAD OF 100 LITRES PER MINUTE?**

The sampling time of operator is reduced of 50% and a double number of tests are performed in comparison with 100 litres per minute.

## **18 - WHY IS SUGGESTED TO ADOPT AN ACTIVE AIR SAMPLER WITH 2 OR 3 ASPIRATING HEADS INSTEAD THAT JUST A SINGLE ONE?**

There are several advantages:

- The use of different media at the same time (e.g.: TSA for bacteria and SDA for moulds)
- The use of the same media for real statistical results
- The use of the same media (e.g.: TSA) at different incubation temperature for bacteria and moulds
- To save operator time
- The application of continuous monitoring during the daily working shift
- To reduce the operator contamination risk because the plates handling is reduced

## **19 - IS IT POSSIBLE TO USE AN ACTIVE MICROBIAL AIR SAMPLER IN EXPLOSION RISK AREAS?**

Yes. The model must be ATEX official certified.

## **20 - WHY ARE 55 MM CONTACT PLATES SOMETIMES ADOPTED INSTEAD OF THE TRADITIONAL 90 MM PETRI DISH FOR MICROBIAL AIR SAMPLING?**


The reason is connected to the fact that some companies prefer to have just a single type of culture plate to reduce the validation of each separated lot.

## **21 - HOW ARE VIRUS COLLECTED FROM AIR?**


The aerosol virus are collected in liquid for rapid analytical identification using the AIR.BIO ONE VIRUS AIR SAMPLER.

# PERFORMANCES

## REDUCE SAMPLING TIME



## GLP-GMP CFR21




**BARCODE IDENTIFICATION 1D-2D  
FOR OPERATOR, LOCATION, PLATE**

## CALIBRATION SYSTEMS



## CASCADE PASSWORDS



**DATA INTEGRITY**



# PERFORMANCES

## ATEX - SHOCKPROOF - IP65



EXPLOSION PROOF  
INSTRUMENTS PRODUCED WITH  
COMPONENTS AND PROCESS  
ATEX CERTIFICATION



TECHNOPLASTIC SHOCKPROOF  
BODY WITH ANTIBACTERIAL  
TREATMENT



WATER & DUST  
PROTECTION IP 65

## INDUCTION - BLUETOOTH



BASE STATION WITH BATTERY CHARGER  
90/260V, AUTOMATIC VOLTAGE SELECTOR,  
CHECK CALIBRATION SYSTEM (OPTIONAL)

**NO PLUG OR  
EXTERNAL CONNECTIONS**

## "AT REST" "IN OPERATION" "END"

UP TO 70 CYCLES  
FOR 1000  
LITRES/HEAD



CONTINUOUS SAMPLING

## FLEXIBILITY

### ASPIRATING HEAD OPTIONS

DISPOSABLE  
TRANSPARENT  
TECHNOPOLYMER  
CERTIFIED STERILE  
"DAILY SHIFT" HEAD



STAINLESS STEEL  
AISI 316  
HEAD



REUSABLE  
TECHNOPOLYMER  
AUTOCCLAVABLE  
HEAD



CONTACT PLATE  
55 MM



PETRI PLATE  
90 MM



FAST CLOSING  
BAYONET






3 DIFFERENT BAYONET  
STYLE ASPIRATING HEADS



SUITABLE FOR  
55 MM CONTACT PLATE OR  
90 MM PETRI PLATE

# FROM SPONTANEOUS GENERATION TO ACTIVE AIR SAMPLER

The Greek philosopher Aristotle (384–322 BC) was one of the earliest scholars who articulated the theory of spontaneous generation, the notion that life can arise from no living matter.



Italian Francesco Redi (1626–1697) and Lazzaro Spallanzani (1729-1799) performed experiments to refute the idea of spontaneous generation.



Luis Pasteur (1822-1895) set experiments that disproved irrefutably the theory of spontaneous generation.

Today the active air sampling confirms the Redi, Spallanzani and Pasteur conclusion theory.







# INDEX

Products	Pag.
TRIO.BAS™ CONTINUOUS AIR MICROBIAL MONITORING	1
TRIO.BAS™ AGC MICROBIAL AIR SAMPLER	7
TRIO.BAS™ MINI	9
TRIO.BAS™ MONO CABLE	11
TRIO.BAS™ MONO FILTER	13
TRIO.BAS™ MONO INDUCTION	15
TRIO.BAS™ AIRBIO ONE CABLE	17
TRIO.BAS™ AIRBIO ONE FILTER	19
TRIO.BAS™ AIRBIO ONE RAPID-VIRUS	21
TRIO.BAS™ DUO CABLE	29
TRIO.BAS™ DUO INDUCTION	31
TRIO.BAS™ AIRBIO DUO CABLE	33
TRIO.BAS™ AIRBIO TRIO CABLE	35
TRIO.BAS™ ATEX (EXPLOSION PROOF)	37
TRIO.BAS™ ISOLATOR RABS	39
TRIO.BAS™ SATELLITE WITH TRI CLAMP CONNECTION	44
TRIO.BAS™ MULTIFLEX 1	45
TRIO.BAS™ MULTIFLEX 1+2	47
ASPIRATING FUNNEL SYSTEM TRI CLAMP	49
ASPIRATING HEAD SYSTEM TRI CLAMP	53
SELFTEST SYSTEM	57
VERITEST SYSTEM	59
TRIO.GAS SYSTEM	61



# INDEX

Products	Pag.
VERIGAS	63
FLUGAS SYSTEM	65
TRIO.CPS	67
TRIO SETTLE	69
DAILY SHIFT HEAD	71
AS SOFTWARE	75
BAS SOFTWARE	77
P.A.C.A.S. SYSTEM	79
CFU PHOTO CAMERA	81
ACCESSORIES	84



# CONTINUOUS AIR MICROBIAL MONITORING IN CLEANROOM

THE MOST COMPLETE RANGE OF MICROBIOLOGICAL AIR SAMPLERS  
FOR SAMPLING UP TO 9 HOURS IN CONTINUOUS

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED




AIRBIO TRIO CM



AIRBIO ONE CM



AIRBIO DUO CM



TRIO.BAS RABS ISOLATOR CM + 1  
SATELLITE CM



TRIO.BAS RABS ISOLATOR CM + 3  
SATELLITES CM



TRIO.BAS RABS ISOLATOR CM + 2  
SATELLITES CM

ORUM INTERNATIONAL introduces a new family of air samplers: AIRBIO and TRIO.BAS RABS ISOLATOR for Continuous Monitoring. The instruments are designed to run at 25 l/m and to be compliant with a continuous environmental monitoring in grade A and B area (GMP Annex 1).

Different solutions available thanks to the great flexibility of the system and the multi-aspirating heads.

It is possible to monitor in continuously from a minimum of 1 hour up to 9 hours (with a 3 aspirating heads air sampler).





## AIRBIO ONE CM

### IDENTIFICATION CODES

Code	AIRBIO ONE for CONTINUOUS MONITORING
------	--------------------------------------

853K AIRBIO ONE CM 25 l/m PETRI with cable PACK (25 litres/min flow rate) (\*)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 1 s/s aspirating head (code 830) with s/s cover head (code 465), 1 cable for transfer data and 1 robustus carrying case.

## AIRBIO DUO CM

### IDENTIFICATION CODES

Code	AIRBIO DUO for CONTINUOUS MONITORING
------	--------------------------------------

854K AIRBIO DUO CM 25 l/m PETRI with cable PACK (25 litres/min flow rate) (\*)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 2 s/s aspirating head (code 830) with s/s cover head (code 465), 1 cable for transfer data and 1 robustus carrying case.



## AIRBIO TRIO CM

### IDENTIFICATION CODES

Code	AIRBIO TRIO for CONTINUOUS MONITORING
------	---------------------------------------

855K AIRBIO TRIO CM 25 l/m PETRI with cable PACK (25 litres/min flow rate) (\*)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 3 s/s aspirating head (code 830) with s/s cover head (code 465), 1 cable for transfer data and 1 robustus carrying case.

# TRIO.BAS RABS ISOLATOR CM AND SATELLITE CM

## IDENTIFICATION CODES

Code	RABS ISOLATOR for CONTINUOUS MONITORING
850K	TRIO.BAS RABS ISOLATOR CM 25 l/m PETRI with 1 SATELLITE PACK (25 litres/min flow rate) (*)
(*) each PACK consists of: 1 TRIO.BAS RABS ISOLATOR with battery charger, 1 calibration certificate, 1 s/s satellite, 1 s/s aspirating head (code 830) with s/s cover head (code 465), 1 cable connection (5 mt) (code 265), 1 cable for data transfer and 1 robustus medium carrying case.	
Code	SATELLITE UNIT for CONTINUOUS MONITORING
851K	CONTINUOUS SATELLITE UNIT Petri PACK
Code	EACH PACK CONSIST OF:
830	S/S ASPI HEAD PETRI CM
465	COVER HEAD s/s to protect ASPI HEAD
265	Cable connection between the Satellite and the command unit - 1 set of 5 meters



**TRIO.BAS RABS ISOLATOR CM + 1  
SATELLITE CM**



**TRIO.BAS RABS ISOLATOR + 2  
SATELLITES CM**



**TRIO.BAS RABS ISOLATOR CM + 3  
SATELLITES CM**

# REMOTE SYSTEM



## REMOTE FUNNEL

### IDENTIFICATION CODES

Code	REMOTE FUNNEL SYSTEM
181	REMOTE FUNNEL SYSTEM in s/s AISI316 with tri clamps connection, funnel and bell - carrying case included



## REMOTE HEAD

### IDENTIFICATION CODES

Code	REMOTE HEAD SYSTEM
192	REMOTE HEAD SYSTEM Petri 90 plate in s/s AISI316 with tri clamps connection - carrying case included

Code	ACCESSORIES for REMOTE SYSTEM
140	s/s valve for REMOTE SYSTEM with 2 TRI clamps connection
184	silicone tube with 2 TRI clamps connection - 50 cm lenght
185	silicone tube with 2 TRI clamps connection - 100 cm lenght
186	silicone gasket for TRI clamp
187	s/s cap clamps
188	s/s 90° elbow short tube with 2 TRI clamps connection - wheelbase 40 mm
193	s/s tri clamps with silicone gasket
194	s/s tube with 2 TRI clamps connection - 100 mm lenght

Code	ACCESSORIES for REMOTE SYSTEM
194/1	s/s tube with 2 TRI clamps connection - 200 mm lenght
194/2	s/s tube with 2 TRI clamps connection - 500 mm lenght
195	s/s 90° elbow long tube with 2 TRI clamps connection - wheelbase 89 mm
197	s/s hexagonal pipe hook with inserts - 100 mm
189	s/s wall connection with 2 TRI clamps connection (wall < 5mm)
198	s/s wall connection with 2 TRI clamps connection - size 80 mm lenght (wall > 5mm)
199	s/s floor pole - base diameter 250 mm - lenght 1 mt



# QUESTIONS AND ANSWERS

## “ENVIRONMENTAL CONTINUOUS MONITORING IN CLEANROOM”

### 1 - WHY IS IT NECESSARY IN GRADE “A” AND “B” OF CLEANROOM TO ADOPT A VIABLE CONTINUOUS MONITORING?

This is what is written in the GMP Annex 1 at the Chapter n.9.

9.24 Continuous viable air monitoring in grade A (e.g. air sampling or settle plates) should be undertaken for the full duration of critical processing, including equipment (aseptic set-up) assembly and critical processing. A similar approach should be considered for grade B cleanrooms based on the risk of impact on the aseptic processing. The monitoring should be performed in such a way that all interventions, transient events and any system deterioration would be captured and any risk caused by interventions of the monitoring operations is avoided.

### 2 - WHICH IS THE RISK OF CONTINUOUS PROLONGED MONITORING WITH “ACTIVE AIR SAMPLING”?

The culture agar medium is progressively dehydrated and it loses the nutrient characteristics for specific nutrients for different micro-organisms. The consequence is that the microorganisms can have difficulty to multiply and the risk is the difficulty to have the correct colonies count.

### 3 - WHAT IS POSSIBLE TO DO FOR THE LIMITATION OF THE EXCESS DEHYDRATION OF THE CULTURE MEDIUM?

It is necessary to adopt an air sampler with a lower air flow rate (e.g.: 25 l/m) and to adopt Petri dishes with at least a volume of 30 ml of culture agar medium.

### 4 - WHY FOR CONTINUOUS MONITORING IS IT SUGGESTED TO ADOPT THE AIR SAMPLER WITH MORE THAN A SINGLE ASPIRATING HEAD?

Several are the reasons.

1. It is possible to divide the working shift time in different periods and in case of microbial contamination, it is possible to identify when the fact happened.
2. The operator handles the air sampler only at the beginning and at the end of activity and therefore the time is reduced with less possibility of contamination.
3. The adoption of the air sampler with satellites gives the possibility to monitor different risk positions.

### 5 - WHICH IS THE SUGGESTED ACTIVE AIR SAMPLER IN ISOLATOR?

The suggested air sampler is the “aspirating funnel” model because the housing of the culture plate is outside from the isolator and therefore the space is saved and the handling with the isolator gloves is eliminated an alternative is the “aspirating head” model.

### 6 - IS IT SUGGESTED TO ADOPT 2 OR 3 ASPIRATING HEADS?

It depends how long is the working shift. Up to four hours the 2 aspirating heads model is o.k. With longer time, the 3 aspirating heads are suggested.

## **7 - WHY TO USE "DAILY SHIFT" STERILE ASPIRATING HEADS?**

The adoption of the sterile "Daily Shift" irradiated certified and triple packed for use in critical cleanroom, RABS and Isolators aspirating heads allows a full traceability throughout the whole supply chain. The certification of sterility simplifies the activity of the laboratory and is convenient for the document inspection.

## **8 - IT IS CONVENIENT TO USE THE "DAILY SHIFT" FROM A COST POINT OF VIEW?**

The calculation of cleaning, packaging for autoclaving, sterilisation process, certificate of sterilisation, storage of documentation, risk of un-correct certification, paper working confirm the advantage of the adoption of "Daily Shift" aspirating head.

They furthermore allow time saving during periods of unusually heavy workloads.

## **9 - WHY IS THE TRANSPARENCY OF "DAILY SHIFT" ASPIRATING HEAD IS AN IMPORTANT ISSUE?**

The transparency of the "Daily Shift" gives the possibility to the operator to control the presence of the culture agar plate.

## **10 - WHICH IS THE PROTOCOL FOR THE VALIDATION OF THE CONTINUOUS ENVIRONMENTAL MONITORING?**

The result of counted cfu of a reference micro-organism on TSA culture plate, with and without the aspiration cycle are compared to validate the same results (with a statistical microbial contest).

## **11 - HOW OFTEN SHOULD THE AIR SAMPLER WITH AN AIR FLOW RATE OF 25 L/M CALIBRATED?**

It is requested the traditional annual control.

## **12 - WHAT ABOUT THE COST OF THE MATERIAL FOR CONTINUOUS MONITORING IN COMPARISON WITH THE MEMBRANE SYSTEM OR OTHER "CLOSED" SYSTEM?**

The final cost of the membrane filter or other "closed" system is higher in comparison with a normal culture agar plate as used with the TRIO.BAS. This is an important fact to consider when you take the decision to invest money for the continuous monitoring air sampler.

## **13 - IS IT POSSIBLE THE VISUAL RECORD OF THE COLONIES (CFU) COUNTED IN THE PETRI DISH?**

Yes, using the photo-camera, in combination with the BAS software, the image of culture plate is stored for evaluation and printing.

# AGC MICROBIAL AIR SAMPLER

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



A COMPLETE INSTRUMENT FOR ALL THE MICROBIAL ENVIRONMENTAL MONITORING NEEDS TEST

**Air**

Sampler

**Gas**

Sampler

**C**

Check Sampler



**AGC**

ACTIVE AIR SAMPLING



GAS SAMPLING



AIR FLOW CHECK SAMPLING



AGC air sampler consists of one instrument to sample air at 100 l/m.

The aluminum bell placed on the aspirating head allows to sample gas at 100 l/m.

The technopolymer yellow bell is used to check the precision of the instrument's flow rate.



## Description

The body of the AGC instrument is designed to stand much more stable on the surface. It is also easier for the end user to read the keyboard and, as the aspirating head is vertical, to place the plates or the bell on it. Primary applications are for pharmaceutical aseptic filling suites, cleanrooms, biotech, IVF clinics, operating theaters, hospital, pharmacies, blood banks.

**It is used for active microbial air sampling.**

**If used for sampling the gas,** the instrument works by measuring the pressure's variation generated by air sampler while air is aspirated through an aluminum bell. A differential pressure sensor measures that variation and compares it with the reference values.

The digital control unit works as a flow meter before the gas passes through the microbial sampler.

**When used to check the status of the calibration,** the instrument runs by measuring the pressure's variation generated by air sampler while air is aspirated through a technopolymer bell placed on the head of the sampler. The differential pressure sensor measures that variation and compares it with the reference values.

At the end of the test, the instrument gives the result of the test: OK if the air sampler is still calibrated, WARNING or ERROR if the air sampler is not calibrated within initial calibration specifics.

- The data can be always transferred to a PC if a dedicated software is installed.
- The transfer data can be made via Bluetooth or via cable.
- The battery is recharged by a power cable connected directly to the air sampler.
- It is possible to work either in manual or automatic mode.
- When used as active air sampler the use of optional sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- It is compact and easy to transfer.
- SOP (Standard Operating Procedure) available from Application Notes.

## Performances

- Technopolymer shockproof body with antibacterial performances of surfaces
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contaminations
- For gas test, the aluminum bell chamber and valve, the regulator are in AISI 316 rated stainless steel. The bell's gasket is in silicon. All parts are autoclavable
- For calibration check the bell chamber is in technopolymer
- Volume of aspirated air: 100 litres/min both for air sampling and gas sampling
- Selected volumes from 30 to 2.000 litres and 17 prefixed programs
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 volt 50/60 hz or by rechargeable battery (inside the air sampler)
- Battery cycles autonomy: 30.000 litres
- Language: English, French, German, Spanish, Italian
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration of 50 users and 50 places
- Delay time, fraction time and fraction number
- Bluetooth connection or cable for data transfer to tablet or PC (with AS SW or BAS SW installed)
- Automatic next calibration reminder
- CE mark
- Continuous/trending analysis according to the USP
- Input pressure of compressed air or gas: 1 ÷ 6 bar
- Suitable for 90 mm Petri dishes or 55 mm Contact plates
- Size instrument: 303x158x135 H mm – weight 1830 gr.
- Stainless steel bell diam. 80x200 h mm – weight 1200 gr.
- Technopolymer bell chamber diam. 100x110 h mm. – weight 300 gr.
- Built in ISO 9001 premises
- IQ, OQ, PQ documentation are available
- Data integrity CFR 21 and GAMP5 (with BAS SW)
- Compliant with ISO 14698-1, EN 17141, ISO 8537-7 and FDA

## Identification Codes

Code	AGC MICROBIAL AIR SAMPLER PACK (*)
670K	AGC MICROBIAL AIR SAMPLER 100 Contact PACK with cable (100 l/m flow rate) for Contact Plate
671K	AGC MICROBIAL AIR SAMPLER 100 Petri PACK with cable (100 l/m) for Petri dishes

(\*) each pack consists of: 1 air sampler with Bluetooth and power battery charger with 1 s/s aspiration head and 1 s/s cover head, 1 compact stainless steel bell chamber with a valve and a regulator for gas pressure, 1 technopolymer bell for calibration check, 1 autoclavable tube for connecting the control unit to the bell, 1 calibration certificate, 1 cable for data transfer, 1 robust carrying case

# TRIO.BAS MINI

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Single aspirating head air sampler with Bluetooth  
and cable for charging



- 100 or 200 litres per minute flow rate model
- Battery charger via cable (110/240 volt)
- Bluetooth for remote control of the air sampler
- Battery cycle autonomy 30.000 litres
- Manual mode only
- AISI 316 rated stainless steel aspirating head with quick bayonet closure and identification number
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Validated according to "EN 17141"

## DESCRIPTION

- Main customers for TRIO.BAS MINI are agro-food industries, dairy, catering, HACCP, beverage, cosmetic, sewage treatment plant, outdoor environment, primary and secondary schools: mainly customers who make few numbers of controls.
- Bluetooth capability allows remote control of the air sampler only. Data transfer via Bluetooth is not an option with the MINI.
- The sampler is free of any external plugs and it is IP65 certified.
- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator and sampling times.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).


## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Battery cycle autonomy: 30.000 litres
- IP65 certified protection from dust and water
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Language: English, French, German, Spanish, Italian
- Operative aspirating cycles: manual
- Memorized data: up to 1.000 samples
- Delayed, remote, start, simultaneous or interval sampling
- Automatic next calibration reminder
- CE mark
- Dimension: 33x16x15h cm
- Weight: 1.430 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	TRIO.BAS MINI PACK (*)
152K	TRIO.BAS MINI 100 Contact PACK (100 litres/min flow rate)
153K	TRIO.BAS MINI 100 Petri PACK (100 litres/min flow rate)
162K	TRIO.BAS MINI 200 Contact PACK (200 litres/min flow rate)
163K	TRIO.BAS MINI 200 Petri PACK (200 litres/min flow rate)

(\*) each PACK consists of: 1 TRIO.BAS MINI with a battery charger, 1 calibration certificate, 1 s/s ASPI head with s/s cover head, 1 light carrying case.



Stainless steel aspirating head with quick bayonet closure



Easy manipulation



Vertical hook to fix the sampler in vertical position on a cart with wheels



# TRIO.BAS MONO cable

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Single aspirating head air sampler with Bluetooth  
and cable for charging



- 100 or 200 litres per minute flow rate model
- Battery charger via cable (110/240 volt)
- Bluetooth for data transfer
- Manual or automatic mode
- Official calibration certificate
- Operator/Administrator cascade passwords
- May be used for compressed gas testing in conjunction with gas system
- Stable on a work surface in a vertical position without the use of any external support
- Validated according to "EN 17141"

## DESCRIPTION

- Main customers are pharmaceutical aseptic filling suites, cleanroom, biotech, IVF clinic, operating theatre, hospital pharmacies, blood banks, clinic, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring and health authorities.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop via a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.
- It is possible to work either in manual or automatic mode.
- The sampler is IP65 certified.
- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile Daily Shift aspirating heads reduces the risk of contamination.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Cycles battery autonomy: 30.000 litres
- IP65 certified protection from dust and water
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Programmable for compressed gases/air
- Data integrity CFR 21 and GAMP5
- CE mark
- Dimension: 33x16x15h cm
- Weight: 1.440 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	TRIO.BAS MONO cable PACK (*)
211K	TRIO.BAS MONO 100 Contact PACK with cable (100 litres/min flow rate)
212K	TRIO.BAS MONO 100 Petri PACK with cable (100 litres/min flow rate)
213K	TRIO.BAS MONO 200 Contact PACK with cable (200 litres/min flow rate)
214K	TRIO.BAS MONO 200 Petri PACK with cable (200 litres/min flow rate)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 1 s/s ASPI head with s/s cover head, 1 cable for data transfer, 1 robustus carrying case.



TRIO.BAS MONO with cable



Stable on a work surface in a vertical position



Easy manipulation

# TRIO.BAS MONO Filter

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Single aspirating head air sampler with Bluetooth,  
cable for charging and HEPA Filter



- 100 litres per minute flow rate model
- Battery charger via cable (110/240 Volt)
- Bluetooth for data transfer
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- HEPA filter
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Use in horizontal or vertical position without any external support
- Validated according to "EN 17141"



## DESCRIPTION

- Main customers are pharmaceuticals, cleanrooms and biotech industries. The HEPA filter allows to capture the particulates.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop. PC or laptop request a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.
- It is possible to work either in manual or automatic mode.
- The sampler is IP65 certified.
- The battery is recharged by a power cable connected directly to the air sampler.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- HEPA filter
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Cycles battery autonomy: 30.000 litres
- IP65 certified protection from dust and water
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 sampler
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Dimension: 33x16x15h cm
- Weight: 1.440 gr
- Built in ISO 9001 premises


## IDENTIFICATION CODES

Code	TRIO.BAS MONO Filter PACK (*)
170K	TRIO.BAS MONO Filter 100 Petri PACK with Cable (100 litres/min flow rate)
171K	TRIO.BAS MONO Filter 100 Contact PACK with Cable (100 litres/min flow rate)
173	Sterile HEPA Filter - highly efficient, single use, bidirectional bacterial/viral removal filter. (10xbox)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 1 box Hepa Filter, 1 s/s ASPI head with s/s cover head, 1 robust carrying case, 1 cable for transfer data.



Easy manipulation



TRIO.BAS MONO Filter detail



# TRIO.BAS MONO induction

TRIO.BAS™

Single aspirating head air sampler with Bluetooth and battery induction charger

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



- 100 or 200 litres per minute flow rate model
- The base station induction charger could be replaced by "SELFTEST SYSTEM" to verify the correct flow rate at regular intervals
- VERITEST is another independent option device to monitor the correct flow rate at regular intervals
- May be used for compressed gas testing in conjunction with gas system
- Manual or automatic mode
- Stable on a work surface in a vertical position without the use of any external support
- AISI 316 rated stainless steel aspirating head with quick bayonet closure and identification number
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Validated according to "EN 17141"

## DESCRIPTION

- Main customers are pharmaceutical aseptic filling suites, cleanroom, biotech, IVF clinic, operating theatre, hospital pharmacies, blood banks, clinic, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring and health authorities.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop via a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler

and a smartphone or tablet (Android version) and then to a PC or laptop.

- It is possible to work either in manual or automatic mode.
- The air sampler is IP65 certified.
- The battery is recharged by a base station induction charger without any cable connection between the air sampler and the charger.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Cycles battery autonomy: 30.000 litres
- IP65 certified protection from dust and water

- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or sampling interval
- Bluetooth connection for data transfer
- Automatic next calibration reminder
- Programmable for compressed gases/air
- Data integrity CFR 21 and GAMP5
- CE mark
- Dimension: 33x16x15h cm
- Weight: 1.440 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	TRIO.BAS MONO induction PACK (*)
200K	TRIO.BAS MONO 100 Contact PACK (100 litres/min flow rate)
201K	TRIO.BAS MONO 100 Petri PACK (100 litres/min flow rate)
205K	TRIO.BAS MONO 200 Contact PACK (200 litres/min flow rate)
206K	TRIO.BAS MONO 200 Petri PACK (200 litres/min flow rate)

(\*) each PACK consists of: 1 air sampler with Bluetooth, 1 calibration certificate, 1 base station induction charger, 1 s/s ASPI head with s/s cover head, 1 robustus carrying case.



Easy manipulation of the aspirating head



SELFTEST SYSTEM: To verify the air sampler correct flow rate at regular interval



TRIO.BAS with induction battery charger

# AIRBIO ONE cable

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



A stationary and portable single aspirating head with Bluetooth capability and cable for charging



- Stable on a work surface in a vertical position without the use of any external support
- Ideally positioned for monitoring laminar air flow systems
- 100 or 200 litres per minute air flow
- Clear easy-to-read visual display
- Battery charger via cable (110/240 volt)
- Bluetooth capability for data transfer
- Cable for data transfer
- Sample up to 30.000 litres on a fully charged battery
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Up to 1.000 memorized data
- Save sampling time (200 litres) by doubling the aspirated volume of air
- May be used for compressed gas testing in conjunction with GAS SYSTEM (optional)
- Validated according to "EN 17141"



## DESCRIPTION

- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation in compliance with quality standards and QM/GMP.
- Primary applications are for pharmaceutical aseptic filling suites, cleanrooms, biotech, IVF clinics, operating theaters, hospital, pharmacies, blood banks, clinical microbiological labs, HVAC systems monitoring, environmental labs.
- An optional barcode module, via the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code).
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop, via a dedicated software installed.
- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.
- The battery is recharged by a power cable connected directly to the air sampler.
- It is possible to work either in manual or automatic mode.
- The 200 lts/min air flow reduces the operators and time sampling.
- The use of optional sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- Possibility to have the version with external HEPA filter that allows to retain the particulates.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according to EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contaminations
- Volume of aspirated air: 100 or 200 lt/m
- Selected volumes from 30 to 2.000 litres and 17 prefixed programs
- The aspirating chamber is suitable for 55 mm Contact plate or 90 mm Petri dish
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 volt 50/60 hz or by rechargeable battery (inserted inside the air sampler)
- Battery cycles autonomy: 30.000 litres
- IP65 protection certificate from dust and water
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder.
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to the USP
- Dimension: 15x15x24h cm
- Weight: 1600 gr.
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	AIRBIO ONE cable PACK (**)
445K	AIRBIO ONE 100 Contact PACK with cable (100 litres/min flow rate)
446K	AIRBIO ONE 100 Petri PACK with cable (100 litres/min flow rate)
447K	AIRBIO ONE 200 Contact PACK with cable (200 litres/min flow rate)
448K	AIRBIO ONE 200 Petri PACK with cable (200 litres/min flow rate)

(\*\*) each pack consists of: 1 air sampler with Bluetooth and power battery charger, 1 s/s aspiration head with s/s cover head, 1 calibrate certificate, 1 cable for data transfer, 1 robustus carrying case.



Easy plate manipulation



# AIRBIO ONE Filter

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



A stationary and portable single aspirating head air sampler with Bluetooth, cable for charging and HEPA Filter



- 100 litres per minute flow rate model
- Battery charger via cable (110/240 Volt)
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- Bluetooth for data transfer
- HEPA filter
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Validated according to "EN 17141"

## DESCRIPTION

- Main customers are pharmaceuticals, cleanrooms and biotech industries. The HEPA filter allows to capture the particulates.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop. PC or laptop request a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.

## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- HEPA filter
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)

- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.
- It is possible to work either in manual or automatic mode.
- The sampler is IP65 certified.
- The battery is recharged by a power cable connected directly to the air sampler.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.

- Cycles battery autonomy: 30.000 litres
- IP65 certified protection from dust and water
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 sampler
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Dimension: 33x16x15h cm
- Weight: 1.440 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	AIRBIO ONE Filter cable PACK (*)
449K	AIRBIO ONE Filter 100 Petri PACK with Cable (100 litres/min flow rate)
450K	AIRBIO ONE Filter 100 Contact PACK with Cable (100 litres/min flow rate)
173	Sterile HEPA Filter - highly efficient, single use, bidirectional bacterial/viral removal filter. (10xbox)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 1 box Hepa Filter, 1 s/s ASPI head with s/s cover head, 1 robustus carrying case, 1 cable for transfer data.



AIRBIO ONE Filter pack +  
STAND UP Holder(code 377)



AIRBIO ONE Filter cable pack

# AIRBIO ONE LIQUID AIR SAMPLER

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

**AIRBIO ONE LIQUID AIR SAMPLER** is an innovative instrument for airborne viable particle sampling. It is specifically designed for total pathogen monitoring of virus, bacteria, molds and yeast.

The instrument has 2 different applications:

- 1) for collecting virus in a liquid sample for subsequent rapid analytical identification (by PCR)
- 2) for traditional impact on agar culture media plate to count colonies (CFU)

## Description

The collection of micro-organisms and virus on culture media requires time for incubation before results. Instead, much faster results are obtained applying the liquid collection method.

This method has great advantages when applied in the following fields:

- Industrial segment (pharma, agro-food, beverage, etc.) for a quick reaction in identifying the contamination
- Military branches for notifying possible biological attacks
- Hospitals for finding the correct pharmacological product and treatment for patients
- Public institutes (school, restaurant, bar, metro, train, municipality buildings, etc.) for disease/pandemic monitoring

AIRBIO ONE LIQUID AIR SAMPLER also allows a quick assessment of the disinfection protocols' efficacy.

The instrument is the result of the European NATO project EUCLID CEPA 13 ("protection of personnel against pathogenic micro-organisms via air sampling and rapid detection and identification").

## Performances of virus sampling

**Principle:** the volume of air is aspirated and mixed in a pre-analytical liquid.

- The collection liquid system is completely sterilizable as produced in stainless steel
- Collection liquid: water, buffer, nutrient broth
- Quantity of collection liquid: 15 ml
- Battery autonomy: 70.000 litres
- 50 users and 50 places
- Dimension of the instrument: 15x20x33h cm
- Weight: 1.850 gr
- Battery: power supply - 218 VDC 60W
- Operating conditions: T° 0-45°C / RH 10% / 60%

## Performances of active air sampling

**Principle:** the volume of air is aspirated on a culture media

- Air flow 200 l/m
- Battery autonomy: 70.000 litres
- 50 users and 50 places
- Dimension of the instrument: 15x15x24h cm
- Weight: 1.600 gr
- Battery: power supply - 218 VDC 60W
- Operating conditions: T° 0-45°C / RH 10% / 60%

## Identification Codes

Code	AIRBIO ONE LIQUID AIR SAMPLER
2448K	AIRBIO ONE LIQUID AIR SAMPLER PACK consisting of: AIRBIO ONE air sampler 200 l/m, 1 s/s Petri aspirating head with s/s cover head, 1 s/s collection liquid system for virus, 2 conical bottom centrifuge tube PP (225 ml), battery charger, 1 calibration certificate and 1 robust carrying case.
	Accessory
324	Conical bottom centrifuge tube PP - 225 ml (8 x box)
413	Self-sealing sterilization bag for s/s collection liquid system for AIRBIO ONE LIQUID AIR SAMPLER (100 x box)
325	Collecting liquid system by PCR for AIRBIO ONE LIQUID AIR SAMPLER
329	VIRUS rack for conical tube PP - size 205x135x102 mm







1. Remove the s/s aspirating head and cover head from AIRBIO ONE



2. Position the collection liquid system into the AIRBIO ONE LIQUID AIR SAMPLER



3. Take out the conical tube from the collection liquid system



4. Fill in the liquid in the tube



5. The instrument is ready



6. On the air sampler, select VIRUS under SPECIAL SAMPLING menu and press START



7. When the sample is finished, transfer the conical tube to the laboratory for analysis



8. The sample is analysed via PCR system

# Sampling Protocol

1

AIR SAMPLING STRATEGY: place, time, frequency according to validated SOP

2

AIR SAMPLER PREPARATION: volume of air and type of liquid according to the analytical laboratory

---

3

AIR SAMPLER PROGRAMMING: air sampler protocol according to user manual and validated SOP

---

4

SAMPLES COLLECTION: after using a sterile conical tube, place a new sterile tube on the collection liquid system for a new test

---

5

SAMPLES: PCR test, qPCR RT, qPCR

---

6

AIR SAMPLER DECONTAMINATION: use 70% isopropyl alcohol

---

7

NEW SAMPLE: the air sampler is ready for a new test

---



# Virus sampling procedures

Sampling air volumes	The suggested volume of air is 2.000 liters. The air sampler has an autonomy of 70.000 litres, but it could be more whist connected to the main power.
Sampling place	The air sampler should be positioned on the direction of the air flow (between door and windows, close to the air conditioning system of HVAC). In hospital, the sampler should be placed near to patient's bed, in public spaces where there is the highest concentration of people.
Sampling Protocol	A specific SOP (Standard Operative Procedure) needs to be prepared.
Collection liquid	The most common used liquid is phosphate buffer or saline buffer. The volume should be 15 ml. In case the aspirating time is longer, it should be necessary to add sterile water to PBS to avoid salt concentration.
Sample Storage	If the analytical sample is not processed right away, the sample needs to be stored at +4°C.
Sample transfer	The sample should be transferred at +4°C unless different indications from analytical laboratory.
Sample processing	Some protocols indicate a concentration's step (by tangential flow filtration).



AIRBIO ONE for active air sampling



Liquid collected by AIRBIO ONE LIQUID AIR SAMPLER





# AIRBIO ONE LIQUID AIR SAMPLER APPLICATIONS



# Statim

# Lab

*division of Orum International*



## **STATIM LAB is the first ISO/EIC 17025:2017 accredited laboratory that carries out calibrations of microbiological air samplers for flow rates of 100 l/m, 180 l/m and 200 l/m.**

STATIM LAB is an Italian metrological laboratory, a division of the ORUM INTERNATIONAL company.

The laboratory was accredited ISO/EIC 17025:2017 by the American company Perry Johnson Laboratory Accreditation, Inc. (PJLA).

It is the first ISO 17025 accredited laboratory in Italy that carries out calibrations of 100/180/200 l/m flow rates of microbiological air samplers.


ISO/IEC 17025 accredited calibration is the most correct option if the microbiological air sampler is used in regulated industries (Pharma, Medical device, ect...) and where quality standards require calibration in accordance with the most stringent regulations.

STATIM LAB is a division of ORUM INTERNATIONAL, manufacturer of the TRIO.BAS air samplers.

In addition to the TRIO.BAS air microbiological samplers, the laboratory is capable of performing professional calibration services for instruments of any brand, whether they are currently in production or out of production.

The experience in the production, calibration and service of microbiological air samplers has allowed the development of a new system – “patent pending” – with very high precision, capable of directly measuring the air flow of instruments and with a reduced measurement uncertainty.

STATIM LAB guarantees calibration periods of 2/3 days.



**STATIM LAB**  
division of Orum International S.r.l.

via Novara 89, 20153 Milano  
<https://www.statimlab.com/>



# TRIO.BAS DUO cable

TRIO.BAS™



Two aspirating heads air sampler with Bluetooth and cable for charging



- 100 or 200 litres per minute flow rate model
- Bluetooth for data transfer
- Cable for battery charger
- Cable for data transfer
- Stable on a work surface in a vertical position without the use of any external support
- AISI 316 rated stainless steel aspirating head with quick bayonet closure and identification number
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Up to 1.000 memorized data
- Use more than one different culture media at the same time
- Saved sampling time by doubling the aspirated volume of air
- Validated according to "EN 17141"

## DESCRIPTION

- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation and in compliance with quality standards and QM/GMP.
- Main customers are pharmaceutical aseptic filling suites, cleanroom, biotech, IVF clinic, operating theatre, hospital pharmacies, blood banks, clinic, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring and health authorities.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader are transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code).
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.

- It is possible to work either in manual or automatic mode.
- The sampler is IP65 certified.
- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- The possibility to use 2 different aspirating heads allows the use of 2 different culture media at the same time or the ability to sample BEFORE (at rest), DURING (in operation) and at the END of each processing cycle.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Battery cycle autonomy: 70.000 litres
- IP65 certified protection from dust and water
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Programmable for compressed gases/air
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to the USP
- Dimension: 26x28x15h cm
- Weight: 1.630 gr
- Built in ISO 9001 premises


## IDENTIFICATION CODES

Code	TRIO.BAS DUO cable PACK (*)
222K	TRIO.BAS DUO 100 Contact PACK with cable (100 litres/min flow rate)
223K	TRIO.BAS DUO 100 Petri PACK with cable (100 litres/min flow rate)
231K	TRIO.BAS DUO 200 Contact PACK with cable (200 litres/min flow rate)
232K	TRIO.BAS DUO 200 Petri PACK with cable (200 litres/min flow rate)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 2 s/s ASPI head with s/s cover head, 1 cable for transfer data, 1 robustus carrying case.



TRIO.BAS DUO with cable



Stainless steel aspirating head with quick bayonet closure



Stable on a work surface in a vertical position



# TRIO.BAS DUO induction

TRIO.BAS™



Two aspirating heads air sampler with Bluetooth and battery induction charger



- 100 or 200 litres per minute flow rate model
- The base station induction charger could be replaced by "SELFTEST SYSTEM" to verify the correct flow rate at regular intervals
- Stable on a work surface in a vertical position without the use of any external support
- AISI 316 rated stainless steel aspirating head with quick bayonet closure and identification number
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Bluetooth capability for data transfer
- Up to 1.000 memorized data
- Use more than one different culture media at the same time
- Saved sampling time by doubling the aspirated volume of air
- Validated according to "EN 17141"

## DESCRIPTION

- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation and comply with the quality standards and QM/GMP.
- Main customers are pharmaceutical aseptic filling suites, cleanroom, biotech, IVF clinic, operating theatre, hospital pharmacies, blood banks, clinic, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring and health authorities.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code).
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- It is possible to work either in manual or automatic mode.
- The sampler is IP65 certified.
- The battery is recharged by a base station induction charger without any cable connection between the air sampler and the charger
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- The possibility to use 2 different aspirating heads allows to have 2 different culture media at the same time or the ability to sample BEFORE (at rest), DURING (in operation) and at the END of each processing cycle.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Light weight, ergonomic and balanced design to facilitate handling with or without gloved hands
- Technopolymer shockproof body with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Battery cycle autonomy: 60.000/70.000 litres
- IP65 protection certified from dust and water
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection for data transfer
- Automatic next calibration reminder
- Programmable for compressed gases/air
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to the USP
- Dimension: 26x28x15h cm
- Weight: 1.630 gr
- Built in ISO 9001 premises


## IDENTIFICATION CODES

Code	TRIO.BAS DUO induction PACK (*)
220K	TRIO.BAS DUO 100 Contact PACK (100 litres/min flow rate)
221K	TRIO.BAS DUO 100 Petri PACK (100 litres/min flow rate)
225K	TRIO.BAS DUO 200 Contact PACK (200 litres/min flow rate)
226K	TRIO.BAS DUO 200 Petri PACK (200 litres/min flow rate)

(\*) each PACK consists of: 1 air sampler with Bluetooth and 1 calibration certificate, 1 base station induction charger, 2 s/s ASPI head with s/s cover head, 1 robustus carrying case.



TRIO.BAS DUO with induction battery charger



Stainless steel aspirating head with quick bayonet closure



SELFTEST SYSTEM: To verify the air sampler correct flow rate at regular interval

# AIRBIO DUO cable

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



A stationary and portable two aspirating heads air sampler with Bluetooth capability and cable for charging



- 100 or 200 litres per minute flow rate model
- Easy display screen reading
- Battery charger via cable (110/240 volt)
- Bluetooth capability for data transfer
- Cable for data transfer
- Battery cycle autonomy up 70.000 litres
- Stable on a work surface in a vertical position without the use of any external support
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Up to 1.000 memorized data
- More than one different culture media at the same time
- Saved sampling time by doubling the aspirated volume of air
- Validated according to "EN 17141"



## DESCRIPTION

- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation and comply with the quality standards and QM/GMP.
- Main customers are pharmaceutical aseptic filling suites, cleanroom, biotech, IVF clinic, operating theatre, hospital pharmacies, blood banks, clinic, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring and health authorities.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop via a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a

smartphone or tablet (Android version) and then to a PC or laptop.

- It is possible to work either in manual or automatic mode.
- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination and provides the certification of sterility requested by regulatory inspectors.
- The possibility to use 2 different aspirating heads allows to have 2 different culture media at the same time or the ability to sample BEFORE (at rest), DURING (in operation) and at the END of each processing cycle.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- Technopolymer body shockproof with antibacterial performances of surfaces
- Complaint according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 volt 50/60 hz or by rechargeable battery (inserted inside the air sampler).
- Cycles battery autonomy: 60.000/70.000 litres
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to the USP
- Dimension: 27x15x25h cm
- Weight: 2.300 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	AIRBIO DUO cable PACK (*)
479K	AIRBIO DUO 100 Contact with cable (100 litres/min flow rate)
480K	AIRBIO DUO 100 Petri with cable (100 litres/min flow rate)
481K	AIRBIO DUO 200 Contact with cable (200 litres/min flow rate)
482K	AIRBIO DUO 200 Petri with cable (200 litres/min flow rate)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 2 s/s ASPI head with s/s cover head, 1 cable for transfer data, 1 robustus carrying case.



Easy plate manipulation



Pack composition



Secure fixing on tripod



# AIRBIO TRIO cable

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

**A STATIONARY AND PORTABLE THREE ASPIRATING HEADS AIR SAMPLER  
WITH BLUETOOTH CAPABILITY AND CABLE FOR CHARGING**



- 100 or 200 litres per minute flow rate model
- Easy display screen reading
- Battery charger via cable (110/240 volt)
- Bluetooth capability for data transfer
- Cable for data transfer
- Battery cycle autonomy up 70.000 litres
- Stable on a work surface in a vertical position without the use of any external support
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Up to 1.000 memorized data, 100 places identification, 100 operators identification
- More than one different culture media at the same time
- Saved sampling time by doubling the air flow rate from 100 to 200 l/m

## DESCRIPTION

- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation and comply with the quality standards and QM/GMP.
- Main customers are pharmaceutical aseptic filling suites, cleanroom, biotech, IVF clinic, operating theatre, hospital pharmacies, blood banks, clinic, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring and health authorities.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop via a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- It is possible to work either in manual or automatic mode.
- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination and provides the certification of sterility requested by regulatory inspectors.
- The possibility to use 2 different aspirating heads allows to have 2 different culture media at the same time or the ability to sample BEFORE (at rest), DURING (in operation) and at the END of each processing cycle.

## PERFORMANCES

- Technopolymer body shockproof with antibacterial performances of surfaces
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 6.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 volt 50/60 hz or by rechargeable battery (inserted inside the air sampler).
- Cycles battery autonomy: 60.000/70.000 litres
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 100
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to the USP
- Dimension: 27x15x25h cm
- Weight: 2.300 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	AIRBIO TRIO cable PACK (*)
487K	AIRBIO TRIO 100 Contact with cable (100 litres/min flow rate)
488K	AIRBIO TRIO 100 Petri with cable (100 litres/min flow rate)
489K	AIRBIO TRIO 200 Contact with cable (200 litres/min flow rate)
490K	AIRBIO TRIO 200 Petri with cable (200 litres/min flow rate)

(\*) each PACK consists of: 1 air sampler with Bluetooth and battery charger, 1 calibration certificate, 3 s/s ASPI head with s/s cover head, 1 cable for transfer data, 1 robustus carrying case.

# TRIO.BAS ATEX

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



**TRIO.BAS MONO and TRIO.BAS DUO for potentially explosive environments (ATEX)**  
**Microbial Air Samplers are used in Zone 2**  
**Explosion Hazard areas (II 3G Ex iC IIC T4 Gc IP55)**



## PERFORMANCES

- Same performances as TRIO.BAS MONO and DUO
- Built in ISO 9001 premises
- Validated according to "EN 17141"
- The TRIO.BAS ATEX microbial air samplers (MONO, DUO) are built with components and production processes compliant with ATEX (explosion proof) certification

## IDENTIFICATION CODES

Code	TRIO.BAS MONO ATEX PACK(**)
207K	TRIO.BAS MONO ATEX (Explosion proof) Air sampler (100 lts/min) CONTACT 55 plate
208K	TRIO.BAS MONO ATEX (Explosion proof) Air sampler (100 lts/min) PETRI 90 plate
209K	TRIO.BAS MONO ATEX (Explosion proof) Air sampler (200 lts/min) CONTACT 55 plate
210K	TRIO.BAS MONO ATEX (Explosion proof) Air sampler (200 lts/min) PETRI 90 plate

(\*\*) each pack consists of: 1 TRIO.BAS MONO ATEX with Bluetooth, 1 calibration certificate, 1 base station induction charger, 1 s/s ASPI HEAD with s/s cover head, 1 robustus carrying case.

Code	TRIO.BAS DUO ATEX PACK(**)
245K	TRIO.BAS DUO ATEX (Explosion proof) Air sampler (100 lts/min) CONTACT 55 plate
246K	TRIO.BAS DUO ATEX (Explosion proof) Air sampler (100 lts/min) PETRI 90 plate
247K	TRIO.BAS DUO ATEX (Explosion proof) Air sampler (200 lts/min) CONTACT 55 plate
248K	TRIO.BAS DUO ATEX (Explosion proof) Air sampler (200 lts/min) PETRI 90 plate

(\*\*) each pack consists of: 1 TRIO.BAS DUO ATEX with Bluetooth, 1 calibration certificate, 1 base station induction charger, 1 s/s ASPI HEAD with s/s cover head, 1 robustus carrying case.



**TRIO.BAS®**  
BIOLOGICAL AIR SAMPLER

## DIRECTIVES ATEX

In some pharmaceutical and food premises the sampling and analytical instruments may operate in production areas and warehouses where are present inflammable material in form of gas, vapour, fog or mist. An explosive atmosphere may be created in presence of such substances and produces an inflammable mix, in case ignition occurs.

The instruments that operate in such environments need to meet the stringent International Directives referred as ATEX (Atmosphere Explosive) directives which are adopted to protect health and safety in the workplace when there are risks due to the presence of potentially explosive atmosphere.

### EXPLOSION RISK IN PHARMACEUTICAL INDUSTRY

Compressed gases and solvents with specific explosive characteristics are processed during industrial production can be released into the air. The same phenomenon can happen in the analytical laboratory in presence of solvents.

A dangerous explosion occurs when the combustible dust comes into contact with an ignition source like electrostatic discharge, mechanical overheating or chemical reactions or heat from the surrounding surfaces.

### EXPLOSION RISK IN FOOD – AGRO AND BEVERAGE INDUSTRY

The dangers of dust explosions in food and beverage plants have been well publicized, but devastating accidents still take place far too often. Almost a quarter of all reported industrial dust explosions occur in the food and beverage industry. There are hazards that need to be controlled, such as those from flammable gas (fuels for ovens), flammable liquids and vapours (spirit based flavourings and cooking/coating oils).



TRIO.GAS in combination with TRIO.BAS MONO ATEX or DUO ATEX Samplers can be used in Zone 2 - Explosion Hazard areas (II 3G Ex iC IIC T4 Gc IP55)



# TRIO.BAS RABS ISOLATOR

TRIO.BAS™

One external command unit fabricated completely in stainless steel connected to 1 up to 3 satellites with Bluetooth capability and cable for charging

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



- 100 or 200 litres per minute flow rate model
- Battery charger via cable (110/240 volt)
- Bluetooth for data transfer
- Cable connection for satellites from 5 to 20 meters
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Cable for data transfer
- Possibility to add additional satellites
- Use more than one different culture media at the same time
- Saving sampling time by doubling the aspirated volume of air
- Validated according to "EN 17141"

## DESCRIPTION

- The RABS ISOLATOR and the satellites are fabricated in AISI 316 rated stainless steel.
- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation in compliance with the quality standards and QM/GMP.
- Main customers are pharmaceuticals, cleanrooms and biotech industries.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader are transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code).
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.
- It is possible to work either in manual or automatic mode.
- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- This air sampler allows monitoring of separated cleanrooms with a single external command unit. The risk of human contamination is reduced, because the satellite units are permanently inside each cleanroom.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- The possibility to use 1/2/3 different aspirating heads allows to have 1/2/3 different culture media at the same time or the ability to sample BEFORE (at rest), DURING (in operation) and at the END of each processing cycle.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## MORE INNOVATIVE AND ESTABLISHED PERFORMANCES

- AISI 316 rated stainless steel (command unit + satellite units)
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Cycles battery autonomy: 60.000/70.000 litres
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or use of cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to the USP
- Dimension: 25x13x18h cm
- Weight: 3.150 gr
- Built in ISO 9001 premises

## MICROBIOLOGICAL MONITORING WITH RABS ISOLATOR

The TRIO.BAS RABS ISOLATOR is an extremely flexible instrument that can be easily adapted to any different type of isolator or RABS.

There are different satellite versions:

- 1. Standard Stationary Satellite (code 260 - 261).** These satellites are fabricated in AISI 316 rated stainless steel. It is possible to use 90 mm Petri dishes or 55 mm contact plates with stainless steel aspirating heads or sterile technopolymer "Daily Shift" aspirating heads. All types of aspirating head to be ordered separately.  
The compact satellite occupies little space inside the isolator.  
Size: diameter 12 cm, height 12 cm, weight 1.170 gr. (without aspirating head).
- 2. Standard Stationary Satellite with HEPA filter (code 262 263).** This satellite has the same features of the standard model.  
It is supplied with an adapter, positioned on one side, to which a HEPA filter is connected for filtering the expelled air. This format is typically used in cleanrooms. The filter's longevity depends on the frequency of the samples use. Replacement is recommended every 3/6 months. If the HEPA filter becomes clogged before this period, the sampler alarm system warns the operator that the airflow is irregular and therefore it is necessary to replace the filter.  
Laterally there is a holder that allows to position the lid of the Petri dish during the sampling phase and to avoid contamination during handling of the plate.  
Size: diameter 12 cm, height 12 cm, weight 1260 gr.



Standard Stationary Satellite with HEPA filter  
(code 262 263)



s/s TRIO.BAS IN-REST (code 180)


CONNECTIONS BETWEEN THE CONTROL UNIT AND THE SATELLITES

The satellites can be connected to the control unit in different ways:

- **Flexible cable (code 265)** with a max extension of 5 mt. This cable is complete with 4-pin male/female connectors. On request, it is possible to supply cables with a length up to about 20 meters. This cable option is available for all satellites.
- **Stainless steel wall connection (code 267)**. Guarantees a hermetic passage through a wall. The flexible cables are not included. Only for standard satellites.



Male/female connectors



s/s wall connection (code 267)



Flexible cable (code 265)

## IDENTIFICATION CODES

Code	TRIO.BAS RABS ISOLATOR with 1 SATELLITE PACK (*)
268K	TRIO.BAS RABS ISOLATOR 100 Contact with 1 Satellite Pack
269K	TRIO.BAS RABS ISOLATOR 100 Petri with 1 Satellite Pack
270K	TRIO.BAS RABS ISOLATOR 200 Contact with 1 Satellite Pack
271K	TRIO.BAS RABS ISOLATOR 200 Petri with 1 Satellite Pack

(\*) each PACK consists of: 1 TRIO.BAS RABS ISOLATOR with battery charger, 1 calibration certificate, 1 s/s satellite, 1 s/s aspirating head with s/s cover head, 1 cable connection (5 mt), 1 cable for data transfer, 1 robustus medium carrying case.

Code	SATELLITE UNIT PACK (*)
260K	SATELLITE UNIT Contact PACK
261K	SATELLITE UNIT Petri PACK

(\*) each PACK consists of: 1 s/s satellite, 1 s/s aspirating head with s/s cover head, 1 cable connection (5 mt).

(\*\*) second or third satellite to be added to basic sampler.



TRIO.BAS RABS ISOLATOR + 1  
STANDARD STATIONARY SATELLITE



TRIO.BAS RABS ISOLATOR + 2  
STANDARD STATIONARY SATELLITE



TRIO.BAS RABS ISOLATOR + 3  
STANDARD STATIONARY SATELLITE



# APPLICATION NOTES CLEANROOM

STANDARD OPERATIVE PROCEDURE – SOP
Air Passive, Air Active, Compressed Gas, Surface Microbial Sampling - SOP
Photo recording of cfu on culture agar plates - SOP
Trend Analysis - Evaluation of Microbial Environmental Monitoring Results - SOP
Contact plate microbiological surface sampler for clean room validation of “cps surface sampler” – SOP
“Daily Shift” microbial surface sampler - SOP
Periodical monitoring of TRIO.BAS air sampler calibration by “Selftest system” - SOP
Passive Air Sampling – Settle Method - SOP
From paper to paper-less - the correct microbial air and surface sampling in cleanroom and controlled environment according iso standard 14698-1 and iso standard 18593 - SOP
Air flow rate of active microbial air sampler checking at regular intervals - SOP
Cleaning procedures for microbial air sampler use in cleanroom - SOP
Contact plate for surface monitoring use - SOP
Disinfection / sterilization of instrument for compressed gas monitoring – SOP
Compressed Gas Microbial Environmental Monitoring According ISO 8573-7 - SOP
Microbial Air Monitoring in Isolator and RABS - SOP
Direct printing of the sampling parameters of the TRIO.BAS with a Bluetooth printer - SOP
Viable Particles Detection with different Compressed Gas operating in Cleanroom - SOP
Data Integrity Good Practice (DIGP)in Microbial Environmental Monitoring - SOP
TRIO.BAS DUO use in Cleanroom - SOP
TRIO.BAS ISOLATOR use in Cleanroom - SOP

IF YOU WANT TO READ THE FOLLOWING APPLICATION NOTES  
YOU ARE INVITED TO VISIT OUR WEBSITE  
AND DOWNLOAD THE LISTED SOPs

**[WWW.TRIOBAS.COM](http://WWW.TRIOBAS.COM)**

OR SUBSCRIBE TO OUR PERIODIC NEWSLETTERS TO RECEIVE  
UPDATES ON THE PUBLICATION OF NEW APPLICATION NOTES.

**[WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER](http://WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER)**

# TRIO.BAS SATELLITE WITH TRI CLAMP CONNECTION


TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

## FEATURES

- No vacuum involved, but a simple electrical cable
- These satellites are fabricated in stainless steel AISI 316.
- The compact satellite occupies little space inside the isolator.
- Fixed location satellites can be easily mounted both horizontally and vertically.
- The TRI CLAMP coupling system makes satellite easy to assemble and disassemble for maintenance and calibration.
- The electric connection cable to the control units remains protected inside the s/s tube.
- It is possible to use 90 mm Petri dishes or 55 mm Contact plates with stainless steel aspirating heads or sterile technopolymer "Daily Shift" aspirating heads.
- All types of aspirating head to be ordered separately.
- Size: diameter 12 cm, height 12 cm, weight 1.170 gr. (without aspirating head).

Code	IDENTIFICATION CODES
320K	s/s SATELLITE TRI CLAMPS CONTACT PACK for Contact plate - with s/s Aspi Head, s/s Cover Head and cable connector
321K	s/s SATELLITE TRI CLAMPS PETRI PACK for Petri dishes - with s/s Aspi Head, s/s Cover Head and cable connector
193	s/s TRI CLAMPS with silicon gasket
194	s/s TUBE with 2 clamps end 10 cm lenght
198	s/s WALL CONNECTION with 2 Clamps end (wall max 5 cm.)
189	s/s WALL CONNECTION with 2 Clamps end (wall over 5 cm.)
265	Cable connection - 1 set of 5 metres.



# TRIO.BAS MULTIFLEX 1

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



**Stainless steel command unit with one aspiration head incorporated and cable for charging**



- 100 or 200 litres per minute flow rate model
- Battery charger via cable (110/240 volt)
- Bluetooth capability for data transfer
- Cable for data transfer
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- Validated according to "EN 17141"

## DESCRIPTION

- The MULTIFLEX 1 air sampler is made in AISI 316 rated stainless steel.
- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation and comply with the quality standards and QM/GMP.
- Main customers are pharmaceuticals, cleanrooms and biotech industries.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader is transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected is transferred via Bluetooth from the air sampler to a PC or laptop via a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- It is possible to work either in manual or automatic mode.
- The data may be transferred via cable, too. This is helpful for all companies that, due to internal policy, are not allowed to use wireless transfer.
- The battery is recharged by a power cable connected directly to the air sampler.
- While under charging, the air sampler can sample.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- Possibility to add one/two satellites to monitor other points in the cleanroom.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- AISI 316 rated stainless steel
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Cycles battery autonomy: 60.000/70.000 litres
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to USP
- Dimension: 25x13x18h cm
- Weight: 4100 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	TRIO.BAS MULTIFLEX 1 PACK (*)
483K	TRIO.BAS MULTIFLEX 1 100 Contact pack with cable (100 litres/min flow rate)
484K	TRIO.BAS MULTIFLEX 1 100 Petri pack with cable (100 litres/min flow rate)
485K	TRIO.BAS MULTIFLEX 1 200 Contact pack with cable (200 litres/min flow rate)
486K	TRIO.BAS MULTIFLEX 1 200 Petri pack with cable (200 litres/min flow rate)

(\*) each PACK consists of: 1 TRIO.BAS MULTIFLEX with battery charger, 1 aspirating head with s/s cover head, 1 robustus carrying case, 1 cable for data transfer, 1 calibration certificate.



Multiflex 1 with cables for optional satellite



Easy manipulation



# TRIO.BAS MULTIFLEX 1+2

TRIO.BAS™



One external command unit completely fabricated in stainless steel with cable for charging, connected to one fixed and two independent satellites



- 100 or 200 litres per minute flow rate model
- Battery charger via cable (110/240 volt)
- Bluetooth for data transfer
- Cable connection for satellites from 5 to 20 meters
- Cable for data transfer
- Suitable for 55 mm Contact plates or 90 mm Petri dishes
- More than one different culture media at the same time
- Saved sampling time by doubling the aspirated volume of air
- Validated according to "EN 17141"

## DESCRIPTION

- The MULTIFLEX 1+2 air sampler and satellites are fabricated in AISI 316 rated stainless steel.
- This air sampler is especially dedicated to customers who make a large number of controls, in different environments, with a large staff rotation in compliance with the quality standards and QM/GMP.
- Main customers are pharmaceuticals, cleanrooms and biotech industries.
- A barcode module, thanks to the use of a scanner (barcode reader) with Bluetooth, automatically records the operator, place and plates used for the sampling. The data collected by the barcode reader are transmitted directly to the instrument. This solution is useful for those who already use culture plates with barcode or 2-D barcode (QR Quick Response Code). The data collected are transferred via Bluetooth from the air sampler to a PC or laptop. PC or laptop request a dedicated software installed.
- The data is transferred via Bluetooth between the air sampler and a smartphone or tablet (Android version) and then to a PC or laptop.
- It is possible to work either in manual or automatic mode.
- The data may be transferred via cable, too. This is helpful for all

companies that, due to internal policy, are not allowed to use wireless transfer.

- The battery is recharged by a power cable connected directly to the air sampler.
- The 200 lts/min air flow reduces the operator time and the time sampling.
- This air sampler allows to monitor separated cleanrooms with a single external command unit. The risk of human contamination is reduced, because the satellite units are permanently inside each cleanroom.
- The use of sterile "Daily Shift" aspirating heads reduces the risk of contamination.
- The possibility to use 2/3 different aspirating heads allows to have 2/3 different culture media at the same time or to make sampling BEFORE (at rest), DURING (in operation) and at the END of each processing cycle.
- The led on the left side of visual display shows the sampling status according to different numbers of flashes (Waiting delay time, Sampling in progress, Sampling pause).

## PERFORMANCES

- AISI 316 rated stainless steel
- Compliant according EN/ISO 14698-1, GMP and GLP
- Stainless steel aspirating head with quick bayonet closure, identification number and stainless steel cover to prevent contamination
- Volume of aspirated air: 100 or 200 l/m
- Selected volumes from 30 to 2.000 litres and 17 preset programs
- The aspirating chamber is suitable for 55 mm Contact plates or 90 mm Petri dishes
- Auto calibration: power/flow electronic real time control
- Power supply system: the instrument can be charged continuously by AC powered source 110/240 Volt 50/60 Hz or by rechargeable batteries (inserted inside the air sampler)
- Cycles battery autonomy: 60.000/70.000 litres
- Language: English, French, German, Spanish, Italian
- Manual and automatic passwords
- Operative aspirating cycles: manual and automatic
- Memorized data: up to 1.000 samples
- Configuration users and places: 50
- Delayed, remote, start, simultaneous or interval sampling
- Bluetooth connection or cable for data transfer
- Automatic next calibration reminder
- Data integrity CFR 21 and GAMP5
- CE mark
- Continuous/trending analysis according to USP
- Dimension: 25x13x18h cm
- Weight: 4100 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	TRIO.BAS MULTIFLEX 1+2 PACK (*)
474K	TRIO.BAS MULTIFLEX 1+2 100 Contact with cable (100 litres/min flow rate)
475K	TRIO.BAS MULTIFLEX 1+2 100 Petri with cable (100 litres/min flow rate)
476K	TRIO.BAS MULTIFLEX 1+2 200 Contact with cable (200 litres/min flow rate)
477K	TRIO.BAS MULTIFLEX 1+2 200 Petri with cable (200 litres/min flow rate)

(\*) each PACK consists of: 1 TRIO.BAS MULTIFLEX 1 PACK, 1 calibration certificate, 2 s/s satellite, 2 s/s aspirating head with s/s cover head, 2 cable connection (5 mt), 1 robustus carrying case.



Connection between command unit and satellite



Command unit with aspirating head and satellites cables



Easy manipulation

# ASPIRATING FUNNEL SYSTEM TRI CLAMP

TRIO.BAS™

**REMOTE STAINLESS STEEL TRI CLAMP ASPIRATING FUNNEL SYSTEM**  
microbial active air sampling

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



This modular system, customized according to the customer's equipment and needs, is particularly useful when there limited space inside ISOLATOR cabinets or RABS.

TRI CLAMPS are one of the most common types of pipe connections in the Biotech and Pharmaceutical industries.

The command unit (TRIO.BAS, AIRBIO, MULTIFLEX) which is generally outside of the cabinet or outside/inside cleanrooms, programs the sampler. The operator must insert the culture media plate into the aspirating chamber of the instrument.

Validated according to "EN 17141"

## DESCRIPTION



This sampling system is completely fabricated in AISI 316 rated stainless steel. The TRI CLAMP connection simplifies all the cleaning/sterilization operations. The risk of contamination associated with operator activities is reduced. The Remote Aspirating Funnel can be manipulated in multiple directions. The Remote Funnel System is fully sterilizable.

## IDENTIFICATION CODES

### TRIO.BAS REMOTE S/S FUNNEL TRI CLAMP

Code	Individual Components:
181	REMOTE FUNNEL SYSTEM – s/s aspirating funnel end clamp with s/s aspirating bell and tube end clamp to connect the system to an Air Sampler
193	s/s TRI CLAMP with silicone gasket
194	s/s tube with 2 clamps end - 10 cm lenght
194/1	s/s tube with 2 clamps end - 20 cm lenght
194/2	s/s tube with 2 clamps end - 50 cm lenght
186	gasket for TRI CLAMP
187	s/s TRI CLAMP closing cap
188	s/s elbow short tube 90° with 2 clamps end - size wheelbase 40 mm
195	s/s elbow 90° tube with 2 clamps end - size wheelbase 88,90 mm, radius 25,40 mm
198	s/s WALL CONNECTION with 2 clamps end (wall >5 mm)
189	s/s WALL CONNECTION with 2 clamps end (wall <5 mm)

(\*) different sizes on request



s/s remote funnel in horizontal position



Detail of the TRI CLAMP connection





THE REMOTE SYSTEM FOR A VIABLE ROUTINE MONITORING USING ASPIRATING FUNNEL

## APPLICATION NOTES CLEANROOM

IF YOU WANT TO READ THE FOLLOWING APPLICATION NOTES  
YOU ARE INVITED TO VISIT OUR WEBSITE  
AND DOWNLOAD THE LISTED SOPs

**[WWW.TRIOBAS.COM](http://WWW.TRIOBAS.COM)**

OR SUBSCRIBE TO OUR PERIODIC NEWSLETTERS TO RECEIVE  
UPDATES ON THE PUBLICATION OF NEW APPLICATION NOTES.

**[WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER](http://WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER)**

# APPLICATION NOTE

By AROUND LAB NEWS

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

## REMOTE MICROBIAL AIR SAMPLING IN RABS AND ISOLATOR

### KEY WORDS

Agar plate, Aspirating head, Contamination, Isokinetic sampling, Remote sampling, Viable routine monitoring

### INTRODUCTION

The REMOTE SAMPLING is necessary for viable routine monitoring because the introduction of a complete air sampler in the working area reduces the space, could increase the risk of contamination during the manipulation and the exhaust of the sampler can disrupt the laminar flow.

It is therefore suggested to adopt a "REMOTE SYSTEM" for the viable routine monitoring.

### THE REMOTE SYSTEM FOR A VIABLE ROUTINE MONITORING

The REMOTE SYSTEM has three different possibilities.

- The aspirating funnel probe system (also called isokinetic probe) or
- The aspirating head system
- The disposable sterile system



ASPIRATING FUNNEL

The agar plate is outside from the cabinet



ASPIRATING HEAD

The agar plate is inside the cabinet

- The aspirating funnel probe (FIGURE 1) has the advantage of space reduction thanks to the limited dimension, the elimination of agar plate introduction in the cabinet avoiding a possible contamination during the manipulation of operator.
- The negative points are the possible partial viable particles loss from the funnel to command unit inside the transfer tubing.
- The aspirating head (FIGURE 2) has the advantage of the agar plate inside the cabinet for an immediate viable particle collection.
- The negative points are the involvement of operator and the more requested space.
- The disposable system (the agar plate is inside a plastic body connected with a plastic tubing to the command unit) has the advantage to be a closed system ready for use.
- The negative point is the cost and the fact the user is bounded to an unique supplier.

# ASPIRATING HEAD SYSTEM TRI CLAMP

TRIO.BAS™

REMOTE STAINLESS STEEL TRI CLAMP ASPIRATING HEAD SYSTEM  
microbial active air sampling

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



This modular system, customized according to the customer's equipment and needs, is particularly useful when there is limited space inside isolator cabinets or RABS.

TRI CLAMPS are one of the most common types of pipe connections in the food, beverage, biotech, and pharmaceutical industries. This type of connection consists of a gasket compressed between two flanges which are in place with a clamp.

Validated according to "EN 17141".

## DESCRIPTION

This sampling system is completely fabricated in AISI 316 rated stainless steel. Microbial sampling is programmed through the control unit (TRIO.BAS, AIRBIO, MULTIFLEX) which is generally located outside the isolator cabinet, RABS or inside/outside the cleanroom.

The operator must insert the plate with the culture media into the separated aspiration chamber. The TRI CLAMP SYSTEM simplifies all the operations of cleaning/sterilization. The risk of contamination associated with operator activities is reduced.


The aspirating chamber can be manipulated in multiple directions, according to the unidirectional airflow. The Remote Head System is fully sterilizable.

## IDENTIFICATION CODES


### TRIO.BAS REMOTE STAINLESS STEEL TRI CLAMP

Code	Individual Components:
191	REMOTE HEAD SYSTEM CONTACT - s/s aspi head chamber with aspirating head and 1 tube clamp end for 55 mm Contact plate
192	REMOTE HEAD SYSTEM PETRI - s/s aspi head chamber with aspirating head and 1 tube clamp end for 90 mm Petri plate
193	s/s tri clamp with silicone gasket
194	s/s tube with 2 clamps end - 10 cm lenght
194/1	s/s tube with 2 clamps end - 20 cm lenght
194/2	s/s tube with 2 clamps end - 50 cm lenght
186	gasket for TRI CLAMP
187	s/s TRI CLAMP closing cap
188	s/s elbow short tube 90° with 2 clamps end - size wheelbase 40 mm
195	s/s elbow 90° tube with 2 clamps end - size wheelbase 88,90 mm, radius 25,40 mm
198	s/s wall connection with 2 clamps end (wall >5 mm)
197	s/s hexagonal pipe hook with gaskets and adaptor
199	s/s floor pole - diameter 25 cm, 1 mt height
184	silicone tube with 2 clamps ends - 50 cm lenght
185	silicone tube with 2 clamps ends - 100 cm lenght
189	s/s WALL CONNECTION with 2 Clamps end (wall <5 mm)

(\*) different sizes on request



Aspirating head in horizontal position



Aspirating head in vertical position and part of the connection made by a silicone tube





THE REMOTE SYSTEM FOR A VIABLE ROUTINE MONITORING USING ASPIRATING HEAD

---

# APPLICATION NOTES CLEANROOM

IF YOU WANT TO READ THE FOLLOWING APPLICATION NOTES  
YOU ARE INVITED TO VISIT OUR WEBSITE  
AND DOWNLOAD THE LISTED SOPs

**[WWW.TRIOBAS.COM](http://WWW.TRIOBAS.COM)**

OR SUBSCRIBE TO OUR PERIODIC NEWSLETTERS TO RECEIVE  
UPDATES ON THE PUBLICATION OF NEW APPLICATION NOTES.

**[WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER](http://WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER)**

# Statim

# Lab

*division of Orum International*



## **STATIM LAB is the first ISO/EIC 17025:2017 accredited laboratory that carries out calibrations of microbiological air samplers for flow rates of 100 l/m, 180 l/m and 200 l/m.**

STATIM LAB is an Italian metrological laboratory, a division of the ORUM INTERNATIONAL company.

The laboratory was accredited ISO/EIC 17025:2017 by the American company Perry Johnson Laboratory Accreditation, Inc. (PJLA).

It is the first ISO 17025 accredited laboratory in Italy that carries out calibrations of 100/180/200 l/m flow rates of microbiological air samplers.


ISO/IEC 17025 accredited calibration is the most correct option if the microbiological air sampler is used in regulated industries (Pharma, Medical device, ect...) and where quality standards require calibration in accordance with the most stringent regulations.

STATIM LAB is a division of ORUM INTERNATIONAL, manufacturer of the TRIO.BAS air samplers.

In addition to the TRIO.BAS air microbiological samplers, the laboratory is capable of performing professional calibration services for instruments of any brand, whether they are currently in production or out of production.

The experience in the production, calibration and service of microbiological air samplers has allowed the development of a new system – “patent pending” – with very high precision, capable of directly measuring the air flow of instruments and with a reduced measurement uncertainty.

STATIM LAB guarantees calibration periods of 2/3 days.



**STATIM LAB**  
division of Orum International S.r.l.

via Novara 89, 20153 Milano  
<https://www.statimlab.com/>

# SELFTEST SYSTEM

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Induction battery charger incorporates a differential pressure device to monitor the deviation of the flow rate compared to the calibration value.



- Having at least one SELFTEST is recommended as part of your best practice to verify the correct volume of aspirated air (ok - warning - error) when sampling with single or multiple instruments
- Technopolymer aspirating bell chamber
- Base station induction battery charger
- Validated according to "EN 17141"

## DESCRIPTION


- The SELFTEST is a system that, instead of the auto calibration already present in the air sampler, checks the precision of the air flow rate. This check is necessary to avoid invalidation of the tests between annual controls for official certification.
- The bell chamber is connected through a tube to the base station induction battery charger.
- Main customers are GLP/GMP for pharmaceutical, cleanroom and biotech industries.
- The base station includes the power supply for charging the TRIO.BAS air sampler batteries.
- SELFTEST is only suitable for air samplers with induction chargers.
- This system works by measuring the depression generated by the air sampler while air is aspirated through a special lid (aspirating bell) applied to the aspirating head. A differential pressure sensor measures that depression and compares it with the set value stored in the TRIO.BAS sampler under test. The result appears at the end of the test, on the LDC of the TRIO.BAS air sampler. Results are displayed as: OK (the air sampler is calibrated), WARNING or ERROR (the air sampler is not within calibration specifics).
- The result is recorded automatically into the air sampler according to data integrity.
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	SELFTEST SYSTEM (*)
351	Base Station induction battery charger with user SELFTEST check calibration system (100 lt/min) for Contact Plate or Petri 90 mm plate
352	Base Station induction battery charger with user SELFTEST check calibration system (200 lt/min) for Contact Plate or Petri 90 mm plate



Aspirating bell of SELFTEST connected to aspirating head of air sampler



TRIO.BAS MONO under control with SELFTEST



SELFTEST system



# VERITEST SYSTEM

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



To check, at regular intervals, the precision level of air flow rate



- Aluminium aspirating bell chamber
- CE mark
- Command data unit
- Validated according to "EN 17141"

## DESCRIPTION

- The VERITEST is a system that, instead of the auto calibration already present in the air sampler, checks the precision of the air flow rate. This check is necessary to avoid invalidation of the tests between annual controls for official certification.
- Main customers are GLP/GMP for pharmaceutical, cleanrooms and biotech industries.
- The VERITEST is a manual system thus the results need to be reported in a document validated by a quality controller.
- VERITEST is suitable for all air samplers and satellites.
- Sizes: 12,5x5h cm.
- This system works by measuring the depression generated by the air sampler while air is aspirated through a special lid (aspirating bell) applied to the aspirating head. A differential pressure sensor measures that depression and compares it with the set value stored in the TRIO.BAS sampler under test. The results appear at the end of the test, on the LCD of the TRIO.BAS air sampler. Results are displayed as: OK (the air sampler is calibrated), or WARNING or ERROR (the air sampler is not calibrated within calibration specifics).
- Built in ISO 9001 premises.

## IDENTIFICATION CODE

Code	VERITEST SYSTEM
353	VERITEST - check calibration system (100 - 200 lt/min) with power supply 6VDC mA

# APPLICATION NOTES BIOCONTAMINATION CONTROL IN AGRO, FOOD, DAIRY, BEVERAGE, CATERING PREMISES

HVAC (Heating Ventilation Air Conditioning) Microbial Hygiene Inspection
Cleanroom Air Monitoring from the Sterile Air Outlets - SOP
Microbiological Air Sampling in Food and Dairy Industries
Cleanrooms and Associated Controlled Environments – Biocontamination Control – Annex D
Evaluation of Bacteria and Moulds from HVAC System in Food/Dairy Industry
Microbial Air Sampling from HVAC Outlet
ASTM Guidelines for Bioaerosol Monitoring
Best Practices for Microbial Air Monitoring in Food, Dairy, Supermarket Premises
EFFL European Union Report Biocontam. Control in Food Cleanroom
Q&A About HVAC in Food, Dairy, Agro, Beverage Premises

The microbial content of air in production premises, warehouse, food dairy and beverage production sites is pretty important. A regular monitoring of these premises has the purpose to reduce the air bioburden in order to obtain products of better quality, less contamination, longer shelf-life, better customer satisfaction.

**IF YOU WANT TO READ THE FOLLOWING APPLICATION NOTES  
YOU ARE INVITED TO VISIT OUR WEBSITE**

**[WWW.TRIOBAS.COM](http://WWW.TRIOBAS.COM)**

**OR SUBSCRIBE TO OUR PERIODIC NEWSLETTERS TO RECEIVE  
UPDATES ON THE PUBLICATION OF NEW APPLICATION NOTES.**

**[WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER](http://WWW.TRIOBAS.COM/SUBSCRIBE-NEWSLETTER)**

# TRIO.GAS SYSTEM

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Tests microbiological quality of compressed air/gas used in cleanrooms.

To be used with ASPI Gas chamber or TRIO.BAS instrument.

TRIOGAS Samplers can be used: Zone 2 - Explosion Hazard areas (II 3G Ex iC IIC T4 Gc IP55)



- Stainless steel calibrated valve and bell chamber
- Official calibration certificate
- Suitable for 90 mm Petri dishes or 55 mm Contact plates
- A stainless steel calibrated valve guarantees 100 litres per minute flow rate
- Completely fabricated in AISI 316 rated stainless steel
- Each component of the TRIO.GAS adaptor instrument is autoclavable with no disassembly required
- No glass valves or meters to crack or break
- Validated according to "EN 17141"

## DESCRIPTION

- The TRIO.GAS guarantees that product contact air is contamination free within sterile or aseptic manufacturing facilities (e.g. cleanroom).
- Main customers are pharmaceuticals, cleanrooms, food and dairy industries.
- The system is according to ISO Standard 8573-7 and ISO 14698-1.
- Before passing through a TRIO.BAS air sampler, the air flow from the compressed supply is regulated by an autoclavable flow meter.
- All the sampling data is transferred via Bluetooth (if using a TRIO.BAS air sampler) to a tablet/smartphone or via Bluetooth to PC by downloaded dedicated software according to GMP and GLP.
- When TRIO.GAS is used in combination with a TRIO.BAS air sampler, the time is regulated by the software of the air sampler.
- When TRIO.GAS is used in combination with ASPI Gas Chamber, the test is manual and the time is regulated by a timer counter.

## PERFORMANCES

- Made in AISI 316 rated stainless steel
- Totally sterilizable via autoclave with no disassembly required
- Volume of aspirated air: 100 litres/min
- Input pressure: 1 ÷ 6 bar
- The aspirating chamber is suitable for 90 mm Petri dishes or 55 mm Contact plates
- Fully compliant according to ISO 8537-7 and EN/ISO 14698-1 FDA
- Dimension: 40x18x25h cm
- Weight: 5290 gr
- Built in ISO 9001 premises

## IDENTIFICATION CODES

Code	TRIO.GAS SYSTEM
600	TRIO.GAS SYSTEM complete of stainless steel electrovalve, gas connection, stainless steel fixing system for air sampler, 1 carrying case and 1 calibration certificate
452	S/S ASPI GAS CHAMBER - with ASPI HEAD CONTACT plate and cover with adapter for TRIO.GAS
453	S/S ASPI GAS CHAMBER - with ASPI HEAD PETRIplate and cover with adapter for TRIO.GAS

Note: Data must be manually recorded and sample must be timed when using ASPI Gas chamber without a TRIO.BAS instrument.



TRIO.GAS in combination with a TRIO.BAS MONO



TRIO.GAS in combination with an ASPI GAS CHAMBER



# VERIGAS

TRIO.BAS™



Two performances in one instrument (EASYGAS + VERITEST):  
compressed gas test and check of the precision level of flow  
rate at regular intervals



VERIGAS is a complementary instrument for TRIO.BAS air samplers.

VERIGAS consists of:

- a digital control unit with tube to connect the bells chambers.
- a stainless steel bell chamber with an on/off valve and a regulator for gas pressure test
- a technopolymer bell chamber to control the air flow rate during the routine activity

VERIGAS can be used with TRIO.BAS, AIRBIO and MULTIFLEX air samplers.

- Validated according to "EN 17141"

## DESCRIPTION

- The system VERIGAS (EASYGAS+VERITEST) works by measuring the pressure's variation generated by air sampler while air is aspirated through a bell-shaped aspirating chamber (in stainless steel for EASYGAS and in technopolymer for VERITEST) applied to the head of the sampler. A differential pressure sensor measures that variation and compares it with the reference values.

1. For EASYGAS the gas or air flow from compressed supply is regulated by the connected digital control unit which works as a flow meter before passing through the microbial sampler.

2. For VERITEST the result appears at the end of the test. Displayed on the control unit as: OK (the air sampler is still calibrated), or WARNING or ERROR (the air sampler is not calibrated within calibration specifics)

- The EASY GAS checks that product contact air is contamination free within sterile or aseptic manufacturing facilities like Cleanroom
- Any air sampler is easily and aseptically connected to the output of compressed gas

- Before entering the aspiration head of a TRIO.BAS air sampler, the air flow from the compressed supply is regulated by a flow valve
- All the sampling data are transferred via Bluetooth or cable (depending on the TRIO.BAS model used) to a PC by a dedicated software according to GMP and GLP
- 1.000 litres of compressed air/gas impact on the 90 mm agar culture Petri dish or Contact plate to collect the microorganisms
- Compliant according to ISO standard 8573-7 and ISO Standard 14698-1
- It can be used in combination with TRIO.BAS, AIRBIO and MULTIFLEX air samplers
- SOP (Standard Operating Procedure) available from Application Notes

## PERFORMANCES

- For EASYGAS, the bell chamber, the valve and the regulator are in AISI 316 rated stainless steel. The bell's gasket is in silicon. All autoclavable
- For VERITEST the bell chamber is in technopolymer
- Volume of aspirated air: 100 litres/min for EASYGAS – 100/200 litres/min for VERITEST
- Input pressure: 1 ÷ 6 bar
- Compact and easy to transfer

- Suitable for 90 mm Petri dishes or 55 mm Contact plates
- Operates with a 1,5V battery (no main power connection)
- Size: command unit 120x80x80 mm
- Stainless steel bell diam. 80x200 h mm – weight 1200 gr.
- Tecnopolymer bell chamber diam. 100x110 h mm. – weight 300 gr.
- Built in ISO 9001
- I.Q., O. Q., P.Q. documentation

## IDENTIFICATION CODE

Code	VERIGAS
599	VERIGAS – for compressed gas test and air flow rate check - with digital control unit, s/s bell chamber for Easygas, technopolymer bell chamber for Veritest, connection tube and robustus carrying case



VERIGAS (EASYGAS version) mounted on AIRBIO air sampler



VERIGAS (VERITEST version) mounted on AIRBIO air sampler

# FLUGAS SYSTEM

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Microbiological sampler to test quality of Air/Nitrogen/CO2  
compressed gas used in Cleanroom




- This device is developed to test the presence of micro-organisms in compressed gas samples (air, nitrogen, CO2, Argon) supplied from tanks and pipes under pressure.
- The 100 lt/min flow rate of FLUGAS System is guaranteed by a flow regulation system controlled by the operator via a calibrated flow meter.
- The system includes ASPIGAS chamber and the results are recorded manually. Optionally, the test can be performed with a TRIO.BAS 100 lt/minutes and the records are automatic.
- Validated according to "EN 17141"

## DESCRIPTION

- FLUGAS System consists of:
  - stainless steel compact base (1);
  - analogic calibrated vertical flow meter (2);
  - stainless steel AISI316 bell chamber (3);
  - on/off valve and a regulator for flow rate (4);
  - ASPIGAS chamber with Aspirating head (5 – 6);
  - digital timer (7);
  - gas input to be sampled (8)
- The system is compliant with ISO Standard 8573-7 and ISO 14698.
- If the FLUGAS System is used in combination with ASPIGAS chamber, the test is manual and the time is calculated by a digital timer.
- If the FLUGAS System is used in combination with a TRIO.BAS air sampler the whole cycle of the test is programmed automatically, controlled and recorded by the software of the device. The sampling data can be transferred via cable or Bluetooth to a PC with a dedicated software according to GMP and GPL.
- All the components of the FLUGAS System are autoclavable (excluding the vertical flow meter and timer).
- Input pressure:  $1 \div 10$  bar:
  - air flow 100/200 lt/min (200 l/m optional)
- Size 21x17x12H cm. (with flow meter 32,5H cm) Weight: 1,38 Kg. (with flow meter 2 kg).

Code	FLUGAS SYSTEM
597	FLUGAS System Petri - with flow meter 100 l/m, ASPIGAS chamber, aspirating head Petri and digital timer
598	FLUGAS System Contact - with flow meter 100 l/m, ASPIGAS chamber, aspirating head Contact and digital timer

Data must be manually recorded and sample must be timed when using chamber without a TRIO.BAS instrument.



Flugas components



Flugas in combination with AIRBIO ONE filter



Flugas with bell mounted on AIRBIO ONE



Flugas with bell mounted on ASPIGAS CHAMBER



# TRIO.CPS

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

Contact Plate Sampler system for surfaces



- Timer: 10 seconds
- Operating temperature: 0-40°C
- Size: diameter 95x45H mm
- Code 289
- Validated according to "EN 17141"

## DETECTION AND ENUMERATION OF BACTERIA, YEASTS AND MOLDS ON SURFACES


- TRIO.CPS is used to evaluate the correct cleaning of all surfaces in contact with food and dairy products to validate the HACCP and training of the staff.
- TRIO.CPS is used in pharma cleanrooms, biotech plants and healthcare facilities to standardize surface monitoring when using traditional irradiated culture media.
- TRIO.CPS is used in combination with the most commercial Contact plate (e.g. RODAC) diameter 55 mm: the standard weight and 10 seconds with visual display guarantee to meet the ISO 18593 for horizontal surface sampling.
- The metal base guarantees easy sterilization/disinfection.
- The CFU are counted at the end of incubation and the results are reported as CFU/cm<sup>2</sup> or CFU/culture plate.
- Ideal for use with standard culture media: Total Bacterial Count (e.g. Trypticase Soy Agar with Lecithin and Polysorbate 80) and Total Yeast and Molds Count (e.g. Sabouraud Glucose Agar with Lecithin and Polysorbate 80).



- 1 Insert the timer in the TRIO.CPS
- 2 Set for 10 seconds



- 3 Slide the contact plate into position
- 4 Remove the lid from the contact plate



- 5 Press the TRIO.CPS system to the surface to be sampled
- 6 Start the 10 seconds timer and lift up the TRIO.CPS at the conclusion of timed sample



- 7 Apply the lid to the contact plate, remove the plate and transfer to microbiology lab
- 8 Report the results as CFU/plate or CFU/cm<sup>2</sup>

## SETTLE PLATE EXPOSURE UNDER UNIDIRECTIONAL AIR FLOW

Articulate stainless steel support for Passive Air Sampling in table and floor model formats



The settle method is applied with the culture plate in a horizontal or inclined position, according to the direction of the unidirectional air flow.

In a horizontal position, air flow may be obstructed creating a slight turbulence over the surface of agar which could deviate the viable particles from the agar's surface.

A slanted position of the Petri dish allows the air to hit the agar surface and pass by without disturbance.

TRIO.SETTLE is an articulate device developed to standardize, test, and enumerate the presence of microorganisms when applying the passive air monitoring, or settling plate, method.

Passive air sampling provides a quantitative analysis of airborne microorganisms deposited onto the surface of an agar plate over a period of exposure.

TRIO.SETTLE can be fixed on a tripod or other floor support devices.

TRIO.SETTLE is fabricated from AISI 316 stainless steel and fully autoclavable.

Validated according to "EN 17141".



## PASSIVE SAMPLING, COMBINED WITH ACTIVE SAMPLING, IS RECOGNIZED BY REGULATORY ORGANIZATIONS AS USEFUL IN ASSESSING THE MICROBIAL QUALITY OF THE AIR ENVIRONMENT.

### THE DIFFERENCE BETWEEN ACTIVE AIR SAMPLING AND PASSIVE AIR SAMPLING

**The Active Method:** a predetermined volume of air is drawn onto an agar plate at a controlled rate of speed.

An active air sampler provides an estimate of the number of microorganisms impacted onto an agar plate, free-floating or carried on particles, within a given size, within a cubic meter of air.

**The Passive Method:** a static agar plate provides an indication of any microorganisms which might settle out of the air due to gravitational effects.

Open agar plates are exposed to the environmental air for a length of time. The number of agar plates placed in the environment, and the exposure time, depends on the risk evaluation. The suggested media is Tryptic Soy Agar (TSA). The microorganism population that settles on the agar plates are counted and evaluated. Passive sampling is simpler and less expensive compared to active sampling, which requires a device. Passive sampling produces an indication of the settling microbial population; active sampling produces a reliable quantification.

**CFU/Exposure time:** WHO World Health Organization. The tables below are based on World Health Organization's guide on monitoring the environment in vaccine manufacturing facilities with a maximum exposure time of 4 hours (90 mm Petri dishes).

CLEANROOM GRADE	PASSIVE SAMPLING EXPOSURE TIME: 4 HOURS
A	CFU <1
B	CFU 5
C	CFU 50
D	CFU 100

### EN 17141:2020 Cleanroom – Biocontamination Control

EXPRESSION OF RESULTS
The number of CFU per plate, per time (e.g. CFU/settle plate/4 H)

## DESCRIPTION

The TRIO.SETTLE consists of:

- AISI316 stainless steel fabrication
- Upper surface disc to accommodate the open culture plate during sampling
- Lower surface disc to accommodate the lid of the culture plate
- Articulate system to obtain different inclinations to the agar surface related to the direction of unidirectional airflow, avoiding laminar flow turbulence
- System to fix the TRIO.SETTLE to a tripod or a base
- Floor base

## PERFORMANCES

The TRIO.SETTLE ideally standardizes the position of the culture plate and reduces the risk of contamination during sampling.

## REFERENCES

Annex C – C.3.3 Settle plate

## IDENTIFICATION CODE

Code	TRIO.SETTLE
367	s/s TRIO.SETTLE table plate stand with timer – size: diam 12 x 20H cm.
368	s/s TRIO.SETTLE floor plate stand with timer – size: base diam. 25 x 110H cm.



TRIO.SETTLE floor plate stand





# DAILY SHIFT HEAD

TRIO.BAS™


ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

Sterile aspirating head for TRIO.BAS air samplers to be used several times during a daily operation shift.



90 mm PETRI DISH (CODE 341)

55 mm CONTACT DISH (CODE 340)



## DESCRIPTION

- The sterile Daily Shift aspirating heads (DSH - sterile Daily Shift Head) avoid the sterilization process necessary for stainless steel aspirating heads.
- The sterilization is proven by an official certificate. This document is requested by regulatory authorities.
- The double irradiated sterile packaging allows the users to always have aspirating heads ready for use.
- The transparency of the sterile Daily Shift is useful for verifying the culture plate is inserted correctly in the aspirating chamber.
- Main customers for sterile Daily Shift heads are: agro-food industries, dairy, catering, HACCP, beverage, cosmetic, sewage treatment plants, outdoor environments, primary and secondary schools, pharmaceuticals, cleanrooms, biotech, hospitals, clinics, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring, health authorities.
- They are suitable for all TRIO.BAS air samplers.
- Shelf life: 6 years from the date of sterilization.

**Stainless steel Aspirating Head(s) of air samplers must be sterilized daily by autoclaving**


## IDENTIFICATION CODES

CODE	STERILE ASPIRATING HEAD
340	Sterile Daily Shift Aspirating Head for Contact plate 55 mm - in double sterile bag
341	Sterile Daily Shift Aspirating Head for Petri 90 mm plate - in double sterile bag

## BENEFIT TO ADOPTING "DAILY SHIFT" ASPIRATING HEAD

- The double irradiated packaging reduces the risk of contamination during manipulation
- A daily certificate of sterilization is requested by regulatory authorities as part of the usual auditing process
- The use of DAILY SHIFT aspirating heads eliminates the cost of in house sterilization and the task of preparing the certificate of sterility
- DAILY SHIFT heads allow time savings during periods of unusually heavy workloads
- DAILY SHIFT heads are semitransparent to see the culture plate inside the aspiration chamber

Stainless steel Aspirating Head(s) must be sterilized daily by autoclaving



# Why adopt the “Daily Shift” Head?

Take a look at the involved processes  
(one daily operative shift in Cleanroom)

STAINLESS STEEL ASPIRATION HEAD  
90mm PETRI DISH - 55mm CONTACT DISH



**Stainless Steel:** Autoclavable, same lifetime of the air sampler, used for official calibration of the instrument and included in the standard kit.

DAILY SHIFT ASPIRATION HEAD  
90mm PETRI DISH - 55mm CONTACT DISH



**Ready to use in the cleanroom:** Double single packaging and certificate of sterilization.

## STAINLESS STEEL HEAD vs READY TO USE HEAD COSTS COMPARISON

A	B	C	D	E
Cleaning	Protection	Sterilization process	Certificate of sterilization	Doc storage

### STAINLESS STEEL HEAD

**DAILY TOTAL COST**  
**A+B+C+D+E=**  
**2\$+1\$+3\$+5\$+4\$=** **\$15\***

*\*based on average cost (calculated in \$) of operator’s labor and time to perform disinfection and documentation*

### READY TO USE HEAD

**DAILY COST**  
**READY TO USE HEAD**  
**+ (E) =**  
**8\$+2\$ =** **\$10\***

*\*based on average cost (calculated in \$) of aspirating and operator labor and time*

STAINLESS STEEL HEAD	DAILY SHIFT HEAD
\$15 x 20 Working Days = \$300	\$10 x 20 Working Days = \$200
\$300 x 10 Working Months = \$3.000	\$200 x 10 Working Months = \$2.000

**CONCLUSION:** “DAILY SHIFT” heads allow to save money (e.g. \$1,000 in one year). It is always convenient to keep them in case of sterilization document inspection.

# DAILY SHIFT HEAD

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

- Irradiated and triple packed for use in critical cleanroom, RABS and ISOLATORS.
- Individually labeled in class 7 cleanroom.
- Allow a full traceability throughout the whole supply chain.

## STERILE “DAILY SHIFT” HEAD PACKAGING





# “AS” SOFTWARE

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED


This software is simple to use and it is suggested when the users mainly want transfer data to a PC for recording or analysis.

This software allows to download sampling data from a TRIO.BAS instrument to a PC:

- no passwords
- possibility to control the air sampler remotely
- possibility to view all samples, to filter the information, to export in “.pdf, .csv and .asd format (encrypted)
- it works with all our air samplers except for MINI

## IDENTIFICATION CODES

CODE	AIRSAMPLER SOFTWARE
295	“AS Software” - transfer data from instrument or smartphone/tablet to PC by Bluetooth or by cable (for cable models)
300	APP Android “ASAPP” - transfer data from instrument to smartphone/tablet and to PC by Bluetooth



AS Software for managing sampling cycle

Active air  
sampling



Data  
transfer



AS Software  
for managing  
air sampling

**AIRSAMPLER Complete**

Index	Device	Start date time	End date time	Volume	Flow	Place	Head	Mode	Result	Dish	Count	Count user	Image
1	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
2	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
3	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
4	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
5	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
6	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
7	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
8	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
9	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
10	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
11	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
12	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
13	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
14	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
15	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
16	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
17	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
18	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
19	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
20	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
21	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
22	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
23	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
24	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
25	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
26	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
27	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
28	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
29	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
30	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
31	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
32	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
33	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
34	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
35	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
36	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
37	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
38	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
39	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
40	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
41	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	
42	425R034216001	2020/01/01 16:22	2020/01/01 16:28	1000	n/a	n/a	1	n/a	ok	0000000000000000	20	LENN (user not certified) (last of 3)	

Signature:

**Single AIRSAMPLER Report**

Device: 425R034216001

Start datetime: 2020/01/01 16:22

End datetime: 2020/01/01 16:28

Calibration expire date:

Volume: 1000

Sampler user: Manual

Place: Manual

Head: 1

Mode: sampling

Result: ok

Dish: 0000000000000000

Incubation time: 6991 d, 22 h, 49 m

Count: 20 (2019/02/22 15:17) incubation time: 6991 d, 22 h, 49 m

Count user: LENN (user not certified) (last of 3)

Image: unknown

Signature: \_\_\_\_\_ Date, time: 2021/02/02 18:19

**CUMULATIVE**

**SINGLE**

Final report

# “BAS” SOFTWARE

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

US FDA 21 CFR annex 11 requires the producers of drugs, medical devices, active ingredients and other types of industries regulated by FDA implement data control systems inside their organisational processes.

The systems which process data in an electronic format must guarantee the security, the traceability and the integrity of data. They need to satisfy some requirements including the access to the systems through individual credentials (personal password), traceability of events through audit trails, no modifiable reporting data and the documentation of the system's specifications.

BAS Software attends all these requirements. This software is used in combinations with TRIO.BAS air samplers and, optional, with the CFU Photo Camera.

- passwords to log in
- 3 different users: root master, administrator and standard user
- possibility to give different permissions to each single user
- possibility to control the device remotely
- possibility to configure the air sampler, setting users and places directly on the software and transfer this data automatically on the air sampler
- possibility to view all samples, to filter the information, to export in different formats: .pdf, .csv, .xml and .bas (encrypted)
- it is also possible to manage different type of sampling: active, passive, surface and gas (special sampling)

## IDENTIFICATION CODES

CODE	AIRSAMPLER SOFTWARE
296	"BAS Software" - transfer data from instrument or smartphone/tablet to PC by Bluetooth or by cable (for cable models)
302	APP Android "BASAPP" - transfer data from instrument to smartphone/tablet and to PC by Bluetooth



BAS SW for managing sampling cycle



AIRBIO air sampler  
Air Sampler Software

Single BAS Report

Device: 116700010004  
Start datetime: 2019/10/30 18:15  
End datetime: 2019/10/30 18:18  
Calibration expire date:  
Volume: 300  
Sampler user: PHX2  
Place: SALA 0003  
Head: 1  
Mode: sampling  
Result: ok  
Dish: 00000000010020w  
Incubation time: 152 d, 0 h, 34 m  
Count: 20 (2020/03/30 19:53) incubation time: 152 d, 0 h, 34 m  
Count user: alex@pgo (Alex Ligu) (last of 2)  
Image: after counting

Note:


Exported by alex@pgo (Alex Ligu) signature Date: 2020/03/30 19:55

Final report with the plate pictures  
before and after manual count



CFU PHOTO CAMERA







# PACAS SYSTEM

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



P.A.C.A.S. for Microbiological EM according to ISO 14698-1


## FULL ENVIRONMENTAL MONITORING SAMPLING SYSTEM

FROM PAPER  
TO PAPERLESS

### Passive Active Compressed Air Surface

Passive air sampling (natural exposure of agar/settling plates), active air sampling (TRIO.BAS or AIRBIO), compressed gas sampling (TRIO.GAS), and surface sampling (TRIO.CPS) can be used in combination with BAS software and CFU Photo Camera.

The BAS software digitally manages all sampling methods (air + compressed gas + surface) in compliance with CFR21 part 11 directives, producing data of the total sampling plan. Contact plate and Petri plate protocols are available, with the greatest advantage achieved using the same plate format for each sample type: Air (active and passive), compressed gas/air, and surface applications.



Passive air Settle plate sampling



Active air sampling



Compressed gas sampling



Surface sampling



Data transfer



CFU PHOTO CAMERA  
(optional)



BAS Software for managing air,  
compressed gas and surface sampling cycles



Final report including plate pictures  
(air, compressed gas and surface)  
before and after manual colony count

# PACAS SYSTEM


TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



PASSIVE AIR + ACTIVE AIR + COMPRESSED AIR/GAS + SURFACE

## DATA TRANSFER FROM THE AIR SAMPLER TO THE LABORATORY





# CFU PHOTO CAMERA

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

Complementary device to BAS software for picture recording of culture plates



Before counting



After counting



As requested by Good Laboratory Practice and Data Integrity 21 CFR part 11, the BAS Software records culture plates' photos immediately before and after the manual CFU count thanks to the CFU Photo Camera (code 337).

The CFU Photo Camera is connected to the PC with the BAS Software installed and the two pictures are automatically recorded and saved in the software. It is possible to export all detailed information of a single sample in a pdf file or print a paper copy.

## DESCRIPTION

- The CFU Photo camera must be used in combination with the BAS Software.
- Simple use: the operator connects the CFU Photo Camera to the PC (with the BAS SW installed), puts the identified culture dish with the counted signed colonies on the designated illuminated area of the CFU Photo Camera, click "Capture" button on the BAS SW to take the picture and store it.
- The pictures of the culture plate are reported at the bottom of the single report where all information of the sample are detailed.
- The report is stored in the BAS SW and then exported in a pdf file for subsequent evaluation.
- Compliant to ISO 7218.
- The system (BAS SW and CFU Photo Camera) is developed according to the requests of the regulatory inspectors.
- The CFU Photo Camera avoids the presence of a second operator, saving additional costs.
- The price of the CFU Photo Camera is affordable and convenient for all microbiological laboratories which apply with GLP and ISO 17025.



## PERFORMANCES

- Light System: LED technology
- Camera resolution: full HD 1080p
- Construction: technopolymer antibacterial treatment
- Suitable for 90 Petri dishes or Contact plate
- Connection to the printer via BAS SW downloaded to PC
- Power voltage: 24/220 V
- Size 17x18x29 H cm.
- Weight: 2,8 kg.
- CE mark
- Made in Italy

## IDENTIFICATION CODE

Code	CFU PHOTO CAMERA
337	CFU PHOTO CAMERA with cable for transfer data and power supply

### Single BAS Report

<b>Device:</b>	116T200618004
<b>Start datetime:</b>	2019/10/30 18:15
<b>End datetime:</b>	2019/10/30 18:18
<b>Calibration expire date:</b>	
<b>Volume:</b>	200
<b>Sampler user:</b>	PINO
<b>Place:</b>	SALA OGGI
<b>Head:</b>	1
<b>Mode:</b>	sampling
<b>Result:</b>	ok
<b>Dish:</b>	000000000019020a
<b>Incubation time:</b>	152 d, 0 h, 34 m
<b>Count:</b>	20 (2020/03/30 19:53) incubation time: 152 d, 0 h, 34 m
<b>Count user:</b>	alexligu [Alex Ligu] (last of 2)
<b>Image:</b>	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <p>after counting</p>  </div> <div style="text-align: center;"> <p>before counting</p>  </div> </div>

**Note:** Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed eiusmod tempor incididunt ut labore et dolore magna aliqua. Ut enim ad minim veniam, quis nostrum exercitationem ullam corporis suscipit laboriosam, nisi ut aliquid ex ea commodi consequatur. Quis aute iure reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint obaecacat cupidatat non proident, sunt in culpa qui officia deserunt mollit anim id est laborum.

Exported by alexligu [Alex Ligu], signature:

Date, time: 2020/03/30 19:56





## PHOTO REGISTRATION BY “BAS SOFTWARE” OF THE CFU COUNTED IN PETRI DISH

Simple use: the operator has just to transfer the identified culture dish with the counted signed colonies on the reading place of the TRIO.BAS CFU photo camera and press a button on the PC where the BAS software is downloaded. According to the GLP, Data Integrity, 21CFR Part 11, the results of all the analytical tests should not be modified and, in such case, it is necessary to show and justify the reason of modification.


These rules are valid for all the analytical cycle, but in the traditional microbiological laboratory the CFU count represents a “uncontrolled and not confirmed data”. The operator, after counting the colonies (CFU) of each culture plate, registers the number and then discharge the plate. There is no trace of count because the plate cannot be stored.

This fact can be solved by a sophisticated and expensive colony counter or more simply adopting a new TRIO.BAS CFU Count photo camera with the BAS software.

### THE NEW PHOTO-CAMERA-PRINTER

- Photo-camera-printer system for image of culture plate recording and storing at the end of incubation and after colony counting.
- The image of the culture plate is presented at the bottom of the single report that summarizes the complete sampling cycle.
- The documentation can be stored in PC where is download the BAS software and transformed in PDF file for subsequent evaluation.
- In compliance with ISO 7218.
- The System has been developed according to the requests of the regulatory inspectors.
- The TRIO.BAS CFU count photo camera is used in combination with the TRIO.BAS Software.

### Single BAS Report

Device:	116T200618004
Start datetime:	2019/10/30 18:15
End datetime:	2019/10/30 18:18
Calibration expire date:	
Volume:	200
Sampler user:	PINO
Place:	SALA OGGI
Head:	1
Mode:	sampling
Result:	ok
Dish:	000000000019020a
Incubation time:	152 d, 0 h, 34 m
Count:	20 (2020/03/30 19:53) incubation time: 152 d, 0 h, 34 m
Count user:	alexligu [Alex Ligu] (last of 2)
Image:	<div><div>after counting</div><div>before counting</div><div></div></div>

**Note:** Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed eiusmod tempor incididunt ut labore et dolore magna aliqua. Ut enim ad minim veniam, quis nostrum exercitationem ullam corporis suscipit laboriosam, nisi ut aliquid ex ea commodi consequatur. Quis aute iure reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint obsecvat cupidatat non proident, sunt in culpa qui officia deserunt mollit anim id est laborum.

Exported by alexligu [Alex Ligu], signature:

Date, time: 2020/03/30 19:56



# MICROBIAL ENVIRONMENTAL MONITORING ACCESSORIES



TRIO.BAS TABLE HOLDER family

Robust platforms which greatly improve the stability of the air sampler when vertically positioned and reduce the risk of fall or damage.

STAND UP HOLDERS:

Completely made in high density impact resistant technopolymer. The shaped base allows the vertical positioning of the air sampler with the use of just one hand. They can be positioned on any work surface or on all tripods and MINI multi holder cart with wheels. Fixed on vertical tripod, the stand allows the air sampler to be independant from the vertical tripod for charging the battery or removing the plates for sampling.



370



376



377

370 STAND UP HOLDER FOR MINI

Size 150x110x90h mm. – weight 390 gr.

376 STAND UP HOLDER FOR MONO, DUO, TRIO

Size 145x105x85h mm. – weight 256 gr.

377 STAND UP HOLDER FOR AIRBIO

Size 165x165x70h mm. - weight 337 gr.



521



521 - VERTICAL HOOK

Size 13x10x25h cm - weight 300 gr.  
Completely fabricated in AISI 316 rated stainless steel. The shaped base allows the vertical positioning of the air sampler with the use of just one hand. It can be positioned on any work surface or can be fixed on all tripods and MAXI multi holder cart with wheels. Fixed on vertical tripod, the stand allows the air sampler to be independant from the vertical tripod for charging the battery or removing the plates for sampling.



530



530 - WALL/TABLE HOLDER

Size 21x13x16h cm - weight 700 gr. Completely fabricated in AISI 316 rated stainless steel. It can be positioned on any work surface or can be fixed onto a wall to keep the air sampler in the same direction of the air coming out from conditioned port.

## TRIO.BAS FLOOR TRIPOD family


The floor tripods allow for positioning the air samplers higher than work surface and orientating them differently from vertical position.

### 523 - STAINLESS STEEL FLOOR TRIPOD

Adjustable height from 150 cm to 200 cm.

Completely made in AISI 316 rated stainless steel.

A ball joint fixes the air sampler and adjust the orientation of the air sampler. Fabricated in AISI 316 rated stainless steel to avoid particle emissions, this tripod is suitable ideally for cleanrooms. The air sampler can be fixed directly on the tripod or alternatively on a stand up holder (optional) fixed on the tripod.



523



521

523



377

## BALOMETER

### 348 - BALOMETER

BALOMETER - it is a tool to be placed on the aspirating head of the air sampler in order to convey the air coming out of the ventilation systems on the ceilings.

It consists in a funnel holder to be applied by pressure on the aspirating head and a polycarbonate funnel.

The BALOMETER comes in two sizes: diameter max 21 cm - diameter max 32 cm



diameter max 21 cm



diameter max 32 cm



# MAXI FLOOR TRIPOD

TRIO.BAS™

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

Transportable and compact, the MAXI TRIPOD is an innovative floor tripod. Made from light weight aluminum, it is sturdy enough to support TRIO.BAS air sampler at great heights



## PERFORMANCES AND SPECIFICATIONS

- Adjustable height from min.125 cm to max 366 cm
- Completely made in anodized aluminum.
- The air sampler can be fixed directly on the tripod or alternatively on a STAND UP holder (optional) fixed on the tripod.
- Compact 4 section stand with 3 risers
- Smooth cushioning of tube movement when the tripod is raised or lowered protects your fingers
- The wide bases opens up to 107cm to keep the tripod steady

## IDENTIFICATION CODE

Code	TRIO.BAS MAXI FLOOR TRIPOD
387	TRIO.BAS MAXI FLOOR TRIPOD
370	STAND UP HOLDER for TRIO.BAS MINI
376	STAND UP HOLDER for TRIO.BAS MONO S, DUO S, TRIO S
377	STAND UP HOLDER for AIRBIO MONO and DUO

# MULTI HOLDER CART ON WHEELS

## MINI model

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



The stand up holder, which is fixed on the top shelf, keeps the air sampler fixed and stable on the cart



Stand up holder  
fixed on the cart

## DESCRIPTION

- The MINI cart is a good tool for cleanrooms, to move easily the air sampler and keep its accessories all together.
- The cart is equipped with 2 shelves.
- Optional: the stand up holder or the vertical hook holder.

## IDENTIFICATION CODES

Code	MULTI HOLDER MINI
371	MINI MULTI HOLDER cart on wheels - stainless steel AISI 314 and 316 - size 545x385x810H mm - weight 6 kg.

# MULTI HOLDER CART ON WHEELS

## MAXI model

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED



Vertical hook fixed on the cart



370  
376  
377

The stand up holder, which is fixed on the vertical support, helps the operator to remove the air sampler for decontamination or charging the batteries. Alternatively the air sampler can be screwed directly on the vertical support.

### DESCRIPTION

- The MAXI cart is a good tool for cleanrooms, to move easily the air sampler and keep its accessories all together.
- The cart is equipped with 2 shelves.
- The vertical support allows the air sampler to reach higher sampling points. The height can be between 147 cm and 218 cm.
- Optional: the stand up holder or the vertical hook holder.

### IDENTIFICATION CODES

Code	MULTI HOLDER MAXI
372	MAXI MULTI HOLDER cart on wheels - stainless steel AISI 314 and 316 - size 545x385x810H (1470-2180) mm - weight 7,700 kg.



# Stainless Steel MULTIFLEX Holder

TRIO.BAS™

For fixing MULTIFLEX Instruments to the TRIO.BAS s/s Tri-pod (code 523)  
or s/s Multi Maxi cart (code 372)

## FEATURES

- Ideal for placement of MULTIFLEX at various heights and cleanroom points that are not easily accessible.
- Simple, stable design: The MULTIFLEX instrument easily attaches to the Holder without fixing screws.
- Fabricated in AISI 316 Steel in AISI 316 Stainless Steel.

Code	IDENTIFICATION CODES
379	s/s MULTIFLEX HOLDER – sizes: 25x14 cm.
484	s/s TRIO.BAS MULTIFLEX Command unit
372	s/s MULTI MAXI CART on wheels – adjustable height from 100 to 210 cm.





# TRIO.BAS IN-REST

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

## Stainless Steel for support culture plate lid and sampling head during sampling phase

Operators, and their activities, represent the greatest risk of contamination inside cleanrooms.

Operators must work under stringent conditions; required to wear special overalls, gloves, facemasks, and eye masks.

For these reasons, it is advisable to create the best possible working conditions.

Air sampling with microbiological samplers is a complex and delicate operation, whether using an instrument with a single aspirating head, or an instrument with more than one aspirating head.



To facilitate operations and avoid errors and risks of contamination during the sampling phase, the TRIO.BAS IN-REST is particularly useful as a temporary support for culture media plates and/or lids, and instrument sampling heads.

## FEATURES

- Built in AISI 316 stainless steel
- Wedge shape inclined at 150 °
- Autoclavable
- Two-sided Support
- Usable for Petri dishes, Contact plates and aspirating heads
- A single hole on one side, and two holes on the other side, facilitate identification of plates and heads
- Size: 90x120x105H mm.
- Code. 180



# GLP\* EASY RACK AISI 316

ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

TRIO.BAS™

Cleanroom Stainless Steel Rack with handle for stacking, storing, and transporting culture plates (Petri or Contact plate). They are used to maintain the distance inside the incubator according the GLP\* (Good Laboratory Practice).



## DETAILS

- They can be refrigerated, incubated and autoclaved.
- Each rack can hold up to 11 Petri dishes, or 15 Contact plates.
- A wide front opening permits easy access but prevents dishes from sliding out.
- The same rack can be adjusted for Petri or Contact plate.
- Avoids all risk of accident during loading and transport to the incubator for staff bio-hazard safety.
- The weight of only 280 gr. make it light and easy to carry.

\*ISO Standard and GPL request that inside the incubator the culture plates are stacked not more than six and at a distance of 2 cm.



- The small size allows to maximize the space inside the incubator.

- Wire design lets you clearly visualize culture dishes.
- Racks are fully autoclavable.

## IDENTIFICATION CODE

Code

GLP\* EASY RACK AISI 316

175

GLP EASY RACK AISI 316

Petri - capacity up to 11 Petri dish 90 mm. size: diam. 110x215h mm.

Contact - capacity up to 15 Contact plate 55mm. size: diam. 110x215h mm.



# CLEANROOM BAG

BAG LINE

BAG LINE - ORUM INTERNATIONAL  
© ALL RIGHTS RESERVED

## With self-adhesive tape

Sterile bags ideal for safely containing and transporting Petri dishes, Contact plates, Swabs, or other objects, outside of the cleanroom:

- Sterility Traceability.
- Triple-wrapped and individually sealed to ensure sterility.
- Beta sterilization certificate included.
- Color is easy to identify in the cleanroom.
- Write-on oversized area to identify sample contents suitable for use with permanent markers.
- Self-adhesive flap closure prevents contamination after samples are placed inside of the bag.
- Included an easy-open tear line which guarantees the sterility of the single bag.
- The elongated, flat bag design facilitates stacking of multiple bags, easy reading of bag contents, and allows plate lids to stay secured inside the bag during handling.



## IDENTIFICATION CODE

Code	CLEANROOM BAG
381	CLEANROOM BAG – plain sampling bag – polyethylene film – sterilized BETA ray – sizes: 40x13 cm. – 15 x bag – 150 x box.

## HOW TO USE CLEANROOM BAG





CARRYING CASE family

Protective cases are suggested to guard air sampler from damage during movement within the production environment and when sent to service centers for calibration.

The ROBUSTUS carrying cases are a very important accessory because they protect the instruments from any damage during the transport phases and especially when they have to be sent to the metrological laboratories for annual calibrations.


The ROBUSTUS carrying cases are resistant to dust, water (IP67), shocks, temperature changes and vibrations during air and truck transport.



Code	CARRYING CASE
392	ROBUSTUS LARGE carrying case for TRIO.BAS MULTIFLEX 1+2/TRIO.BAS MULTISTATION with 1-3 Satellites - large 56x43x31h cm - medium 56x43x21h cm
394	ROBUSTUS MEDIUM carrying case for TRIO.BAS MULTIFLEX 1+2/RABS ISOLATOR with 1-3 Satellites - 56x43x31h cm
395	ROBUSTUS STANDARD carrying case - 48x38x17h cm
398	ROBUSTUS LARGE carrying case TRIO.GAS - 56x43x31h cm
401	ROBUSTUS MEDIUM carrying case for ARBIO - 56x43x22h cm
403	LIGHT carrying case for VERITEST - 34x30x16h cm
524	LIGHT STANDARD carrying case for TRIO.BAS MINI - 43x35x19h cm

BLUETOOTH PRINTER

DIRECT PRINTING



520 - BLUETOOTH PRINTER

Ultra-light and compact portable Bluetooth printer with high printing autonomy. It has an end paper sensor and prints most popular barcodes. A practical belt hook, a battery charger, one paper roll and USB cable are included.  
Size 11x9x5h cm - weight 450 gr.

421 - ROLL REFILLS

For BLUETOOTH PRINTER - size 57 mm - 10 x box

423 - THERMAL ROLL PAPER

For BLUETOOTH PRINTER - size 57 mm - 10 x box FEATURES - no erasable thermal paper for TRIO. PRINTER - The paper is printed thermally but it keeps the ink during the time

NEW THERMAL PAPER

S.N.:	0000001	S.N.:	0000001
Flow Rate(l/min):	200	Flow Rate(l/min):	200
Date:	15/12/20 12:31	Date:	15/12/20 12:31
Type:	Sample	Type:	Sample
Result:	OK	Result:	Interpreted
User:	Admin	User:	Admin
Place:	400	Place:	400
Volume(1):	50	Volume(1):	50
Sample-1:	200944	Sample-1:	200958
Sample-2:	200948	Sample-2:	200958
Calib.expiry:	11/11/21	Calib.expiry:	11/11/21
Signature		Signature	

This does not fade

This fade




## BARCODE READER



294

291

292



BARCODE READING  
OF PETRI DISH

### 294 - BARCODE READER BLUETOOTH 1D 2D USE

This miniature barcode reader, frequently used in microbiological air monitoring procedures, can help save time, better control the activity of the operators and achieve complete traceability of the sampler tests - size 6x3,5x1,5h cm - weight 45 gr.

### 291 - LOCATION PRESET BARCODE TAG

Size 8,5x5,5 cm - 10 x box

### 292 - USER PRESET BARCODE TAG

Size 8,5x5,5 cm. - 10 x box

## Portable Command Unit (PCU)



301

### 301 - PCU - PORTABLE COMMAND UNIT

301 - PCU PORTABLE COMMAND UNIT. It has a 7" LCD display. The Bluetooth connection with all TRIO.BAS instruments allows download of sampling data. The PCU can also be used to remotely switch on, switch off and pause the air sampler. The data downloaded from the instruments through the PCU can be transferred to a PC where software is installed (i.e. by USB connection). The Portable Command Unit is an ideal instrument to simplify and to facilitate the activity of operator.

BAS SOFTWARE for TRIO.BAS Microbial Air Samplers

295



AS SOFTWARE FOR TRIO.BAS MICROBIAL AIR SAMPLERS

“AS” Software, which stands for Air Sampler Software, is a software that allows the transfer of sampling data from the air sampler to the PC for recording or further analysis. The data transfer is possible via Bluetooth for all air samplers or via USB cable for cable models only.

Code	AIRSAMPLER SOFTWARE
295	“AS” Software - transfer data from instrument to PC via Bluetooth or via cable (for cable models only)
300	APP Android “ASAPP” - transfer data from instrument to Android smartphone / table

296



BAS SOFTWARE FOR TRIO.BAS MICROBIAL AIR SAMPLERS

“BAS” Software is a software which is complaint with CFR 21 Data Integrity as data attends the data standards of ALCOA:

- **A** ttributable as all individuals need to log in with a password
- **L** egible as data is easy to understand and to be recorded permanently
- **C** ontemporaneous as data is recorded when observed
- **O** riginal as data is accessible in the original form
- **A** ccurate as free of errors and complaint with the protocol

The data transfer from air sampler to PC is possible via Bluetooth for all air samplers or via USB cable for cable models only. BAS SW simplifies the entry of user profiles and sampling programs.

Code	AIRSAMPLER SOFTWARE
296	“BAS” Software - transfer data from instrument to PC via Bluetooth or via cable (for cable models only)
302	APP Android “BASAPP” - transfer data from instrument to Android smartphone / table

IQ,OQ, PQ documents



500... AND MORE - IQ, OQ, PQ DOCUMENTS

For industries involved in pharmaceutical and healthcare products or laboratories, equipment quality is very important and even small inconsistencies can generate disastrous results. Installation Qualification (IQ), Operational Qualification (OQ) and Performance Qualification (PQ) are essential components of quality assurance. IQ OQ PQ protocols establish that the equipment, which is installed and used, offers a high quality assurance, so that manufacturing processes will consistently produce products that meet predetermined quality requirements.



# SUPPORT, SERVICE, CALIBRATION AND REPAIRS

ORUM INTERNATIONAL offers an accurate service of maintenance and calibration of the whole range of TRIO.BAS samplers, directly or through the qualified technical assistance centers of its distributors.

In particular, the periodic calibration of the instruments is carried out using best practice metrology methods by highly qualified technicians.

ORUM INTERNATIONAL, as the inventor of portable microbiological samplers with more than 40 years of experience, is able to provide technical support to your staff (QC) to facilitate the implementation of the entire range of TRIO.BAS samplers in your daily testing routine.



We recommend official calibration of the air samplers every 6-12 months. We also recommend a recalibration when the air sampler is potentially damaged, the flow rate is compromised, or any time the firmware is upgraded.

During the recalibration, the air sampler's flow rate is checked to guarantee the value of aspirated air is consistent and the instrument works correctly. A detailed certificate of calibration is subsequently issued after the calibration.

